A guide for using

# Breeze™ RMGI

Resin Modified Glass Ionomer Cement

## PENTRON™

Manufactured By: Pentron Clinical  
1717 West Collins Avenue, Orange, CA 92867 USA  
800-551-0283 / 714-516-7557 / www.pentron.com

©Pentron Clinical 2016. All rights reserved.

Rx ONLY

DMC2753\_Rev.00

### ENGLISH

#### Breeze™ RMGI

Resin Modified Glass Ionomer Cement

#### DESCRIPTION

Breeze™ RMGI is a dual-cure paste/paste resin-modified glass ionomer cement combining the benefits of a resin cement and a glass ionomer cement. It is provided in a dual-barrel syringe delivery system fitted with an automixing tip for ease of application. This cement achieves high bond strength to tooth structure without the need of an etchant or an adhesive. Excess cement can be easily removed by allowing the cement to reach its gel state either through self-curing or tack-curing. Its enhanced radiopacity allows easy identification in x-rays.

#### INDICATIONS FOR USE:

Cementation of:

- Metal-based inlays, onlays, crowns and bridges
- Resin-based inlays, onlays, crowns and bridges
- Ceramic inlays
- Zirconia based all ceramic crowns and bridges
- Fiber, metal, and ceramic posts

#### CONTRAINDICATIONS:

- Pulp capping
- Sensitivity, in rare cases in patients with extreme sensitivity

#### WARNINGS:

Uncured methacrylate resin may cause contact dermatitis and injure the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact.

#### PRECAUTIONS:

1. Hold the restoration with moderate pressure while cleaning excess cement.
2. Do not treat or clean the zirconia-based or metal-based restorations with phosphoric acid or hydrofluoric acid (HF) which can negatively impact adhesion.
3. Do not use provisional cements containing eugenol, as these materials may interfere with the curing of Breeze™ RMGI.
4. Do not open foil package until initial use. Use the product within 4 months of opening the pouch. Always ensure that it is not past the expiration date printed on the syringe label.

#### ADVERSE REACTIONS:

None Known.

#### INSTRUCTIONS:

1. Prepare the tooth. If a liner is needed, use calcium hydroxide or another suitable pulp liner material.
2. Remove the provisional restoration. Clean the preparation with pumice and water. Rinse thoroughly with water after cleaning preparation.
3. Air dry the preparation so that the surface has a somewhat glossy appearance. Do not over dry or desiccate the preparation.
4. Try in the restoration to ensure proper fit and occlusion. Adjust if needed. Thoroughly clean the bonding surfaces of the restoration.
5. Remove the cap from the syringe.

6. If the level of material is not even, bleed excess material from the syringe until the base and catalyst both flow evenly.
7. Attach the mixing tip to the syringe.
8. Dispensing of the material can be done as follows:
  - a. Inlay/Onlay – Dispense the cement directly into the cavity preparation covering all surfaces. Seat the restoration onto the tooth preparation, allowing the cement to flow from all sides.
  - b. Crown – Dispense the cement directly onto the restoration, covering all surfaces. Seat the restoration onto the tooth preparation, allowing the cement to flow from all margins.
  - c. Post – Prepare the post space. Size and fit the post. Dispense the cement onto the post or directly into the canal space, seat the post, and vibrate slightly to avoid trapped air. Allow the cement to slowly flow from the canal space.
9. Once the restoration is seated, remove the excess cement. Excess cement is best removed in its gel state with a scaler. Gel state can be achieved by tack-curing excess with a light for approximately 2-3 seconds, or by allowing the cement to self-cure until it feels rubbery.
10. Material should completely set within 4 minutes and 30 seconds after placement in the mouth.
11. Finish the restoration and check occlusion when material has completely set.
12. Remove the used mixing tip and replace it with the original cap; do not leave the used mixing tip on the syringe.

#### Important Note:

Verify that the laboratory has pre-treated the final restoration in accordance with manufacturer's instructions. Automix tips are for single patient use only, to prevent cross-contamination between patients.

#### Working Characteristics:

Mode of Cure:	Dual Cure
Method of Mixing:	Auto-Mix
Mixing Ratio (base:catalyst)	1:1
Mixing Time	Auto-Mix
Working Time (minimum)*	1 minute 30 seconds
Intra-oral Setting Time (maximum, from initial mix)*	4 minutes 30 seconds

\* These times may vary based on storage conditions, temperature, humidity, and age of product.

#### LIMITED WARRANTY – LIMITATION OF PENTRON'S LIABILITY

Pentron's technical advice, whether verbal or in writing, is designed to assist dentists in using Pentron's product. The dentist assumes all risk and liability for damages arising out of the improper use of Pentron's product. In the event of a defect in material or workmanship, Pentron will replace the defective product or part thereof, or reimburse the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Pentron.

A guide for using

# Breeze™ RMGI

Resin Modified Glass Ionomer Cement

## PENTRON™

Manufactured By: Pentron Clinical  
1717 West Collins Avenue, Orange, CA 92867 USA  
800-551-0283 / 714-516-7557 / www.pentron.com

©Pentron Clinical 2016. All rights reserved.

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