ParaPost®XT™

Threaded Endodontic Post System

Instructions for Use

CE
0344

 COLTENE

EN	Instructions For Use • One-Office-Visit Technique	3
DE	Gebrauchsanleitung • Stiftaufbau in einer Sitzung	6
FR	Mode d'emploi • Technique de reconstitution directe	9
ES	Instrucciones de uso • Técnica en una sola visita.	12
IT	Istruzioni per l'uso • Tecnica diretta unica visita	15
NL	Gebruiksaanwijzing • Directe opbouw techniek	18
SV	Användarinstruktioner • Direktteknik vid ett besök.	21
DA	Brugsanvisning • Eet Kliniksbesøgteknikken	24
FI	Käyttöohjeet • Yhden hoitokerran tekniikka	27
NO	Bruksanvisning • Teknikk for engangsvisitt	30
PT	Instruções para uso • Técnica de uma Consulta.	32
EL	Οδηγίες Χρήσης • Τεχνική μίας επίσκεψης.	35

One-Office-Visit Technique

Intended Use: Metal posts are intended to be cemented into the root canal of a tooth to stabilize and support a restoration.

 CAUTION!	Federal law (United States) restricts this device to sale by or on the order of a healthcare practitioner.
--	--

 NON STERILE	Posts are delivered non-sterile and are to be sterilized before use.
---	--

 WARNING!	Failure to properly follow the directions in this IFU may result in patient injury or compromise the integrity of the restoration.
--	--

 WARNING!	Device should not be used in patients with known material sensitivity to titanium.
--	--

We urge the use of Hygenic® rubber dam with this procedure.

Note: Posts are single use devices and should never be salvaged for re-use. Structural integrity of the post could be compromised

including microscopic damage which could lead to post fracture.

For illustrations see page 40-41.

* Note that titanium posts have virtually the same radiopacity as composite and gutta percha. Careful scrutiny of the radiograph will allow for the differentiation of post from gutta percha and composite.

1. After endodontic therapy has been completed, prepare tooth for restoration as if pulp were intact. The preparation should include at least 1.5 mm of sound tooth structure all around the circumference of the preparation, apical to the core, for desired ferrule effect.* If that amount of tooth structure is not available, it should be created by periodontal crown lengthening and/or orthodontic extrusion.*
2. Use a radiograph to determine the appropriate diameter and depth of post space preparation. Keep in mind that sufficient root wall thickness must be maintained (1 mm minimum)* to prevent perforation or weakening of root walls, and at least 4–5 mm of gutta-percha must be left intact to protect the apical seal.
3. Remove gutta-percha to preplanned depth with a Gates-Glidden drill, Peeso reamer and/or hot instrument. Radiographic verification is recommended. If the canal has been obturated with a silver cone, remove it and reseat

with gutta-percha before preparing the post space.

4. To begin paralleling the post space, select the ParaPostXT Drill, the diameter of which corresponds to the last Gates-Glidden drill used in Step 3 (Fig. 1). (See Fig. A. “Approximate Diameter Comparisons” chart.) To help maintain the preplanned depth of the preparation, ParaPostXT Drills have been pre-marked at 7, 9, and 11 mm from the apical end of the drill.

ParaPostXT Drills can be used either manually with the Universal Hand Driver or with a slow-speed contra-angle (750–1,000 rpm). When a contra-angle is used, the drill must be kept in continuous clockwise rotation until it has been completely removed from the tooth. This will minimize the risk of the drill seizing in the post space. Remove all weakened or unsupported tooth structure. Remove any tooth debris from canal by irrigating post space with water spray.

5. Sequentially step up to the next larger ParaPostXT Drill until the preplanned diameter and depth are achieved (Fig. 2).
6. Select the ParaPostXT post that corresponds to the last drill used to prepare the post space. Before seating the post, it is necessary to determine that the length of the post corresponds to the depth of the newly created post

space. Reinsert the last drill used into the post space to confirm the depth of the preparation, using either the markings on the drill or a rubber stop. Place the post next to the drill so that the base of the post head is slightly above the final depth measurement identified on the drill (Fig. 3). Using a cutting disc, remove that portion of the post which extends beyond the apical end of the drill. Re chamfer the end of the post to its original shape. This should result in a post length which is 0.5 mm shorter than the actual post space (Fig. 4).

7. Align the undercuts on the post head with the prongs on the wrench and place the post head firmly into the hand wrench. Use wrench marked “A” with post sizes 3, 4, and 4.5. Use wrench marked “B” with post sizes 5, 5.5, and 6. For safety, tie dental floss around the cleat on the top of the wrench and secure extra-orally.
8. Insert the post into the post space and slowly engage the threads into the dentin by turning in a clockwise direction (Fig. 5). As resistance is encountered, remove the post by turning the post counterclockwise. Clean and rinse the post space. Clean the post. Reinsert post and continue to thread the post into the dentin until the shoulder gently contacts the tooth (Fig. 5 and 6).
9. Evaluate the post head in relation to the occlusion.

Remove the post by turning in a counterclockwise direction (Fig. 7). Adjust the post head for occlusal clearance (Fig. 8). Clean, rinse and dry the canal with Hygenic Paper Points (Fig. 9). Clean and dry post.

10. Lightly coat the walls of the post space with cement using a lentulo spiral (Fig. 10).
11. Lightly coat the post shaft with cement and insert the post into the post space. Using the hand wrench, slowly engage the threads into the dentin by rotating the post in a clockwise direction. As resistance is encountered, back off a quarter turn and then continue as before (Fig. 11). When the shoulder contacts the tooth, the post is fully seated and should be backed off a quarter turn to further relieve any insertional stress.
12. If necessary, to obtain additional occlusal clearance after the cement has set, make any adjustments needed to the post head while holding the post with a hemostat (Fig. 12).
13. After the cement has set, remove any excess before completing the core and fabricating the final restoration (Fig. 13).

Cleaning and Sterilization

Drills and Metal Wrenches: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush; thoroughly

rinse and dry. Place drills, with a sterilization indicator, into an autoclave bag and seal. Sterilize in autoclave at 132°C for 12 minutes. Sterilizing bags should not touch the walls of the autoclave. Once the sterilization monitor indicates that sterilization was achieved, place autoclave bag, with the drills inside, in storage until ready for use. After use, to decontaminate, soak in a 70% isopropyl alcohol solution for 10 minutes.

Posts: Prior to each use, place posts, with a sterilization indicator, into an autoclave bag and seal. Sterilize in an autoclave at 130° C for 15 minutes, following standard procedures. Alternatively, sterilize in a pre-vacuum chamber at 132° C for 4 minutes, 134° C for 3 minutes or 134° C for 6 minutes.

Universal Hand Driver: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush. Sterilize by autoclave or chemical vapor sterilization.

Warning: DO NOT BEND POSTS

* References and/or literature available upon request.

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©2018
Coltène/Whaledent Inc.
Coltène/Whaledent AG
Coltène/Whaledent GmbH + Co. KG