



# COLTENE

October - December 2024

## Promotions

NEWSLETTER



SIGN UP!

LAUNCH  
PROMOTION

### BRILLIANT Lumina

TOOTH BRIGHTENING



LAUNCH OFFER  
**BRILLIANT Lumina**  
P. 2



FREE EXTRA HEAD  
**CanalPro X-Move**  
P. 12



NOVEMBER OFFER  
**HyFlex EDM - OGSF Sequence**  
P. 16



Professional Tooth Brightening System

# BRILLIANT Lumina

**NEW**

Natural bright smile

**H<sub>2</sub>O<sub>2</sub>**  
was yesterday!

**PAP**  
is the future of  
tooth brightening.

BRILLIANT Lumina offers you an innovative in-office solution for tooth brightening – the first formula without hydrogen peroxide and carbamide peroxide.

This is because BRILLIANT Lumina is based on the active ingredient PAP (phthalimidoperoxy caproic acid). What makes it so special: the unique formula enables the PAP concentration to be kept sufficiently high to achieve a naturally brighter smile during treatment.

No  
hydrogen  
peroxide



LAUNCH  
OFFER

**1+1**  
ON BASIC KITS  
*While stocks last*

ONLY AVAILABLE FROM  
YOUR TERRITORY MANAGER

**BRILLIANT Lumina Basic Kit**  
REF 65002864  
4x Gel, 1x Activator, 1x Patient brochure

LEARN MORE



## BRILLIANT Lumina – an asset for your practice

Satisfied patients – no tooth sensitivity during or after treatment.

Gentle on the tooth structure – analyses of tooth surfaces show that BRILLIANT Lumina does not have an erosive effect on the tooth structure.

Simple treatment protocol

– the uncomplicated application of BRILLIANT Lumina can also be performed by dental practice staff.

### Convincing results

Well-documented patient cases and consistent results confirm the reliability of BRILLIANT Lumina. And patient satisfaction is also impressive.

Before and after images of tooth brightening with BRILLIANT Lumina



*The clinical outcomes depend on the individual circumstances of the patient.*



"I've never had tooth bleaching before because I was worried about unnaturally white teeth and pain sensitivity. However, the gentle method of BRILLIANT Lumina quickly convinced me. The natural brightening of my teeth is exactly what I was looking for – and all without pain during or after treatment."

**Andreas, 27 years**



Light-Curing Single-Component Universal Bond

# ONE COAT 7 UNIVERSAL

THE NEW  
SUPERKRAFT

Strong individually, unbeatable together:  
The **SoloCem** + **ONE COAT 7 UNIVERSAL**  
bonding system

**300**  
Applications  
per bottle

- Self Etch, Selective Etch and Total Etch, a bond for all adhesive techniques
- With activator for chemically curing products
- Exceptional bonding values on enamel and dentine
- Excellent marginal seal and outstanding marginal leakage tightness
- Ergonomically shaped, triangular bottle with dropper for fine dosing

Contains MDP - One bond for all indications



**ONE COAT 7  
UNIVERSAL Refill**  
REF 60019539



**UNIVERSAL ACTIVATOR 3 ml**  
for self and dual cure materials  
REF 7054  
1 x 3 ml Activator,  
1 x 50 pcs. Microblush, black



LEARN MORE



**ONE COAT 7 UNIVERSAL Intro Kit**  
REF 60019538  
1 x 5 ml ONE COAT 7 UNIVERSAL  
1 x 2.5 ml Etchant Gel S  
1 x 20 pcs. Application needles  
1 x 50 pcs. Brushes black  
1 x Mixing well with cover





A partner you can count on – always

# SoloCem

## A host of different indications, yet always the same solution

You believe in speed and efficiency? Then choose SoloCem, because self-adhesive application is not only fast but also safe. However, depending on the clinical case, a stronger adhesive bond may be necessary. Even then, SoloCem is the right choice – use ONE COAT 7 UNIVERSAL as "adhesion booster" in addition.

## Less administration. More time for important things

As you only need one cement for both luting options (adhesive and self-adhesive) with SoloCem, you not only save space in your inventory, but also time spent for the administration of additional products. That means: More time for what is most important – the patient

LEARN MORE



SAMPLES  
AVAILABLE FROM  
YOUR TERRITORY  
MANAGER

For improved bond strength: use with **ONE COAT 7 UNIVERSAL** on enamel and dentine

### SoloCem Refill Dentin

REF 60014062  
1 × 5 ml / 10.5 g SoloCem Dentin,  
5 × Mixing Tip Brown Short Super Fine  
5 × Mixing Tip Brown Short Fine

### SoloCem Refill White Opaque

REF 60014061  
1 × 5 ml / 10.5 g SoloCem White Opaque / as above

### SoloCem Refill Translucent

REF 60019840  
1 × 5 ml / 10.5 g SoloCem Translucent / as above

### SoloCem Intro Kit

REF 60013986  
1 × SoloCem Dentin, 5 ml,  
1 × SoloCem Trans, 5 ml  
10 × Mixing Tip Brown Short Super Fine,  
10 × Mixing Tip Brown Short Fine

Universal submicron hybrid composite

# BRILLIANT EverGlow

## Enduring gloss – made brilliant



96%  
global mark  
obtained  
by publishers  
Dental

With the universal composite BRILLIANT EverGlow highly aesthetic and long-lasting restorations can be achieved effortlessly.

- Exceptional polishability and long-lasting brilliance
- Aesthetic single-shade restorations
- Versatile shade system with three translucency levels
- Good wettability with minimal stickiness to the instrument
- Excellent sculptability and smooth consistency



TERRITORY MANAGER OFFER  
OCTOBER ONLY



# BRILLIANT EverGlow Flow

The ideal complement

- Duo Shade system
- Ideal for cavity linings and small fillings
- Enduring gloss
- Good adaptation properties due to its low viscosity
- High stability and perfect thixotropic
- Absolute level of precision



TERRITORY MANAGER OFFER  
OCTOBER ONLY



### HIGHLY AESTHETIC RESULTS



Initial situation:  
discoloration



Opaque shade  
applied



Final  
restoration





### EverGlow Refills

1 x Syringe, 3g  
OR 20 x Tip, 0,2

	<b>Bleach</b> Syringe > 60019697 Tips > 60019706		<b>Translucent</b> Syringe > 60019704 Tips > 60019718
	<b>A1/B1</b> Syringe > 60019698 Tips > 60019707		<b>Bleach Translucent</b> Syringe > 60019705 Tips > 60019719
	<b>A2/B2</b> Syringe > 60019699 Tips > 60019708		<b>Opaque Bleach</b> Syringe > 60019743 Tips > 60019746
	<b>A3/D3</b> Syringe > 60019700 Tips > 60019709		<b>Opaque A1</b> Syringe > 60019744 Tips > 60019747
	<b>A3.5/B3</b> Syringe > 60019701 Tips > 60019710		<b>Opaque A2</b> Syringe > 60023819 Tips > 60023823
	<b>C2/C3</b> Syringe > 60019702 Tips > 60019716		<b>Opaque A3</b> Syringe > 60019745 Tips > 60019748
	<b>A4/C4</b> Syringe > 60019703 Tips > 60019717		<b>Opaque A4</b> Syringe > 60023820 Tips > 60023824

### EverGlow EcoPacks

<b>A1/B1, Tips EcoPack,</b>	60 x 0.2 g 60024160
<b>A2/B2, Tips EcoPack,</b>	60 x 0.2 g 60024161
<b>A3/D3, Tips EcoPack,</b>	60 x 0.2 g 60024162
<b>A3.5/B3, Tips EcoPack,</b>	60 x 0.2 g 60024163
<b>A1/B1, Syringes EcoPack,</b>	5 x 3 g 60024166
<b>A2/B2, Syringes EcoPack,</b>	5 x 3 g 60024167
<b>A3/D3, Syringes EcoPack,</b>	5 x 3 g 60024168
<b>A3.5/B3, Syringes EcoPack,</b>	5 x 3 g 60024169



LEARN MORE



### EverGlow Flow

1 x 2g Syringe,  
5 x Standard application needles  
5 x Fine application needles 0.4 mm

	<b>Bleach BL</b> Syringe > 60019752
	<b>A1/B1</b> Syringe > 60019753
	<b>A2/B2</b> Syringe > 60019754
	<b>A3/D3</b> Syringe > 60019755
	<b>A3.5/B3</b> Syringe > 60019756
	<b>A4/C4</b> Syringe > 60019757
	<b>Translucent Trans</b> Syringe > 60019758
	<b>Opaque A2/(OA2)</b> Syringe > 60019751
	<b>Opaque A3/(OA3)</b> Syringe > 60023825

# Dual Curing Composite Based Post Cement and Core Build Up Material

LEARN MORE



## ParaCore

3 indications – 1 material

- 3 indications – 1 material
- For post cementation, core construction and cementation of crowns, bridges, inlays and onlays



### ParaCore 5 ml Intro Kit

REF 5885

- 1 x 5 ml ParaCore Dentin,
- 1 x 5 ml ParaCore White
- 1 x 5 ml ParaCore Dentin SLOW,
- 1 x 5 ml ParaCore Trans SLOW
- 3 x 3 ml ParaBond Non-Rinse, Conditioner + Adhesive A & B
- 20 x Mixing Tip Short Super Fine




REF 5853	Dentin Refill,	1 x 25 ml
REF 5854	White Refill,	1 x 25 ml
REF 5886	Automix White Refill,	2 x 5 ml
REF 5887	Automix Dentin Refill,	2 x 5 ml
REF 60011391	Slow Dentin Refill,	2 x 5 ml
REF 60011392	Slow Trans Refill,	2 x 5 ml

## Temporary Restorations

## DuoTEMP

Double security keeps you smiling



- Dual-curing
- High marginal seal due to slight expansion
- Good retention, easy removal
- Contains zinc oxide and fluoride for additional protection
- Eugenol-free - no effect on the curing of composites
- Radiopaque

### DuoTEMP Single Pack

REF 5831 1 syringe of 5 g

### DuoTEMP Eco Pack

REF 7090 5 syringes of 5 g



Debridment / Root canal shaping

# MicroMega HandFiles



- Pre-sterilised, ready-to-use files
- Guaranteed sterilisation in compliance with regulatory standards

FULL RANGE OF FILE SIZES AVAILABLE  
PLEASE ASK YOUR TERRITORY MANAGER

## MicroMega K-Reamers

Pre-sterilised instruments  
Blisters of 6 HandFiles



- REF 60029655 06
- REF 60029658 08
- REF 60029661 10
- REF 60029664 15
- REF 60029667 20
- REF 60029670 25

## MicroMega K-Files

Pre-sterilised instruments  
Blisters of 6 HandFiles



- REF 60029718 06
- REF 60029721 08
- REF 60029724 10
- REF 60029727 15

## MicroMega H-Files

Pre-sterilised instruments  
Blisters of 6 HandFiles



- REF 60029781 08
- REF 60029784 10
- REF 60029787 15

LEARN MORE



# CanalPro SyringeFill

## Smart colour-coded syringes for safe root canal irrigation

Modular system for safe and efficient filling of syringes with endo irrigating solutions. The colour coding of the CanalPro irrigation solution system simplifies handling and ensures safe treatment workflows, even in stressful situations. CanalPro SyringeFill is available in two versions. Once as CanalPro SyringeFill Station with a free-standing station consisting of base and bottle carrier in a stainless steel construction. In addition, it is also available as CanalPro SyringeFill Holder with a bottle carrier for wall mounting in a stainless steel construction.

- Protects against contamination
- Improved safety by colour coding to minimise the risk of syringe mix-up
- Time-saving one-hand use
- Prevents loss of material
- Available as a free-standing station or as a holder with bottle carrier for wall mounting



SyringeFill Station  
free standing - REF 60031347

SyringeFill holder  
wall mount - REF 60031348



# CanalPro Irrigation Solutions

Irrigation solutions are used to help dissolve vital and necrotic tissue, eradicate bacteria and other microorganisms, remove the smear layer, and lubricate the canal system.

As no single irrigation solutions exists that can accomplish all of the above goals, multiple solutions should be used during clinical treatment. Sodium hypochlorite (NaOCl), ethylenediaminetetraacetic acid (EDTA), and 2% Chlorhexidine (CHX) are commonly used for chemical debridement during the process of canal preparation.



**SyringeFill Cap**  
REF 60019662

## CanalPro NaOCl

3% or 6%

For irrigation/debridement of root canals during and after treatment

Available in 3% and 6% formula

500ml - 3% - REF 60019655

500ml - 6% - REF 60019656

## CanalPro EDTA 17%

17% EDTA solution (pH 8.5)

Removes smear layer and dentin mud

Opens dentin tubules for:  
- Disinfection solutions to work better  
- Better adhesion of sealers and obturation materials

500ml - REF 60019651

## CanalPro CHX 2%

2% chlorhexidine digluconate solution

Substantive

Very effective against *E. faecalis* and fungal

500ml - REF 60019658



**NEW**

# CanalPro X-Move

Includes integrated Apex Locator

Already programmed for the new OGSF sequence file

- Continuous rotation and a reciprocation motion for all common file systems
- Wireless endomotor for easy and flexible operation
- Integrated apex locator for safe and efficient treatments
- Insulated contra angle with super mini-head and slim neck provides convenient usage and free field of view
- Intuitive control display for easy handling
- Pre-programmed settings for COLTENE file systems saves preparation time




CanalPro X-Move contra angle super mini-head & slim neck



Normal contra angle head & neck

LEARN MORE



## CanalPro X-Move Motor REF 65002739

Contents:

- 1 × Motor Handpiece
- 1 × Charging Base
- 1 × Contra Angle Head (REF 65002744)
- 1 × Spray Nozzle
- 1 × Measuring Wire
- 4 × File Holder
- 2 × Lip Clip
- 2 × Touch Probe
- 2 × "O" ring
- 1 × Power Adapter
- 1 × Battery



**OFFER**

BUY AN X-MOVE MOTOR AND GET AN

**EXTRA CONTRA ANGLE HEAD FREE!\***

Promo code: REF 60036546

Ask your Territory Manager for a demonstration



Digital assistance system for canal preparation

# CanalPro Jeni

Autonomous "navigation" in the endo canal owing to revolutionary software

The new CanalPro Jeni endomotor has a patent-pending digital assistance system for canal preparation, which uses complex algorithms to control the file movement at millisecond intervals. The rotary motion and speed are continuously controlled based on current intensity and torque.

Different routes, obstacles and situations CanalPro Jeni is named after its developer Prof. Dr. Eugenio Pedullà and guides you safely through the traffic like a navigation system. As such, the Jeni system automatically adapts to the individual root canal anatomy, working in unison with the mechanical and chemical preparation every step of the way.

LEARN MORE



## CanalPro Jeni, complete set REF 60023659

includes: 1 × Control unit, 1 × Motor with apex measuring contact, 1 × Fully insulated angled handpiece 1:1, 1 × Bluetooth foot switch, 1 × Power supply, 1 × Apex cable set, 1 × Holding bracket for apex cable

MicroMega

# Dual Move

## Dual Move Motor

REF 52002004

Option: 1 × Contra-angle head  
REF 51200158

LEARN MORE



- Cordless motor, well-balanced in hand and ambidextrous
- Miniature head, 340° rotative contra-angle
- Continuous rotation or reciprocation modes
- Fully customisable
- Accurate apex localisation
- Apex location – combined use

# Electrical Discharge Machining

## HyFlex EDM

Stays on track

Manufactured from Nickel-Titanium alloy with Controlled Memory (CM) heat treatment, HyFlex EDM files allow multiple regeneration and reuse. Due to their innovative manufacturing process using electric discharge machining HyFlex EDM files have a unique surface and stand out by their durability and high fracture resistance.

- Treats 80% of all cases with the one file
- Designed for successful centered canal preparation
- Multiple regeneration and reuse
- Specially hardened for safer application
- High flexibility
- Extreme resistance to fracture

***"I have been amazed at how predictably I can achieve excellent outcomes"***

*"Since my conversion to COLTENE HyFlex EDM files, I have been amazed at how predictably I can achieve excellent outcomes, even when taking on challenging teeth. I have gained confidence with the files' ability to effortlessly navigate sclerosed and curved canals. It is a simple system that can be used in almost all cases and it gives you a feeling of reassurance when you have such a high-quality file."*

*Their fracture resistance makes it possible for me to treat complex cases in my general practice setting, which is great for my patients.*

*Additionally, I have found the system to be very efficient, as it reduces the need for additional files, making endodontics more affordable."*



**Dr. Orry Torsney**  
Principal of Central Park Dental,  
Mullingar, Co Westmeath.

LEARN MORE



**FULL RANGE OF FILE SIZES AVAILABLE**  
PLEASE ASK YOUR TERRITORY MANAGER



**OFFER**

**2+1**

ON REFILLS

MINIMUM ORDER 4

TERRITORY MANAGER OFFER  
NOVEMBER ONLY

- Pre-Sterile**  
for safer application
- High Fracture resistance**  
Uniquely hardened surface reduces risk of file separation
- Variable Cross Section Design**  
Contributes to high breakage resistance
- Triangular cross section at top**
- Cutting Efficiency**  
Hardened surface plus controlled memory improves cutting efficiency
- Trapezoidal cross section in middle**
- Highly Flexible**  
Controlled memory allows high flexibility and pre-bending
- Quadratic cross section at tip**



Evo Endo and COLTENE Academy bring you

# Everything Endo

## THE A-Z OF ENDODONTICS IN PRACTICE

An intensive 3 day hands-on course  
with Dr Dhiraj Arora

**Day 1 Friday 8th November 2024**

**Day 2 Saturday 9th November 2024**

**Day 3 Saturday 23rd November 2024**

9am - 4.30pm

**etc venue, 50-52 Chancery Lane, London WC2A 1HL**  
(Nearest Tube Stations: Chancery Lane or Holborn)



### Dr Dhiraj Arora

BDS MJDF RCS(Eng) MSc (Endo) PG Cert CE

Dr Dhiraj Arora graduated from Guy's, King's and St Thomas' Dental Institute, completing his Bachelor of Dental Surgery in 2007. In 2009, he went on to complete his postgraduate examinations and gain The Diploma of Membership of the Joint Dental Faculties at The Royal College of Surgeons of England.

Following his passion in endodontics, Dhiraj completed his Masters in Endodontic Practice at Barts and The Royal London. His postgraduate training has equipped him with the latest knowledge and skills required to provide high quality and evidence based endodontic treatment in practice. He is a member of the British Endodontic Society and the European Society of Endodontology. Dhiraj's primary research interest is in the use of CBCT and digital imaging in Endodontics.

The aspects he enjoys most about his job are helping patients and educating other dental practitioners. His keen interest in teaching, led him to gain a Postgraduate Certificate in Dental Education from Oxford Brookes University. He has been actively involved in the mentorship of newly qualified dentists and is an experienced educational supervisor with over 10 years' experience within the London Deanery. This, along with being an Honorary Lecturer in Endodontics (QMUL) allows him to combine his enthusiasm for Endodontic treatment and teaching.

An intensive 3-day, hands-on course, designed to help the general dental practitioner learn and understand modern techniques required for daily endodontic procedures, using the latest materials and equipment

We spread the course over 3 days so that you can consolidate your learning in between sessions. Also you get much more from three separate dedicated days and you don't need to find three consecutive days to be away from work.

### OBJECTIVES

Following this course delegates should be able to:

- Provide an accurate endodontic diagnosis and understand treatment planning considerations for assessing the restorability of a tooth and endodontic treatment
- Gain confidence and feel comfortable with rubber dam isolation
- Achieve predictable endodontic access and canal orifice location, appreciating the anatomy of a tooth
- Understand and implement modern mechanical preparation protocols using multiple systems
- Carry out effective and safe endodontic irrigation
- Obturate canals in a predictable fashion, free of voids and to an adequate length
- Restore an endodontically treated tooth with a fundamental understanding of when deep marginal elevation is an effective treatment modality

**Cost:** £1,795 inc. vat

**GDC Outcomes:** C and D


**For course information,  
availability and to book**





**NEW**

# HyFlex EDM – OGSF Sequence







## Would you like fast, safe and predictable? Discover the new sequence solution!

The new HyFlex EDM OGSF sequence includes a new Orifice Opener and Glidepath file for opening the glide path, as well as a Shaping file for shaping the root canal over the full length. The defined treatment procedure makes the preparation easy and keeps the learning curve short. Due to the optimized match of tip sizes and tapers, designed in cooperation with Prof. Pedullà, the files provide a smooth feeling when shifting from one file to another and help you to keep control over the treatment.

### Features and Benefits

- Efficient treatment with high cutting performance
- High fracture resistance and flexibility
- Centered preparation of the canal
- Only one shaping file
- Sterile blister
- 2 lengths: 21 and 25 mm

### Shaping

HyFlex EDM OGSF Set, 4 pcs	21 mm	25 mm
OGSF Sequence	60036141	60036142
HyFlex EDM OGSF Refills, 3 pcs	18 mm	21 mm
Orifice Opener 18/.11	60036140	-
Glidepath File 15/.03	-	60023493
Shaping File 18/.045	-	60036136
Finishing File 30/.04	-	60036138
HyFlex EDM Finishing File, 3 pcs	21 mm	25 mm
Finishing File 40/.04	60022271	60019597
Finishing File 50/.03	60022272	60036598
Finishing File 60/.02	60022274	60036599

### Drying and Obturation

Size / Taper	HyFlex EDM Paper Points, 100 pcs	HyFlex EDM Guttapercha Points, 60 pcs
30/.04	60036503	60036504
40/.04	60022676	60022708
50/.03	60022741	60022709
60/.02	60022740	60022710



TERRITORY MANAGER OFFER  
NOVEMBER ONLY

LEARN MORE



# ParaPost System

A complete range of posts for direct and indirect indications

## ParaPost System Refills

REF PF1613	ParaPost Fiber White Posts Size 3	5 pcs
REF PF16145	ParaPost Fiber White Posts Size 4.5	5 pcs
REF PF1615	ParaPost Fiber White Posts Size 5	5 pcs
REF PF16155	ParaPost Fiber White Posts Size 5.5	5 pcs
REF PF1616	ParaPost Fiber White Posts Size 6	5 pcs
REF PF1713	ParaPost Fiber Lux Posts Size 3	5 pcs
REF PF1714	ParaPost Fiber Lux Posts Size 4	5 pcs
REF PF17145	ParaPost Fiber Lux Posts Size 4.5	5 pcs
REF PF1715	ParaPost Fiber Lux Posts Size 5	5 pcs
REF PF17155	ParaPost Fiber Lux Posts Size 5.5	5 pcs
REF PF1716	ParaPost Fiber Lux Posts Size 6	5 pcs
REF PF18145	ParaPost TaperLux Posts Size 4.5	5 pcs
REF PF1815	ParaPost TaperLux Posts Size 5	5 pcs
REF PF18155	ParaPost TaperLux Posts Size 5.5	5 pcs
REF PF1816	ParaPost TaperLux Posts Size 6	5 pcs
REF P6830	ParaPost XT Titanium Posts (.036"/0.90mm)	10 pcs
REF P6840	ParaPost XT Titanium Posts (.040"/1.00mm)	10 pcs
REF P6845	ParaPost XT Titanium Posts (.045"/1.14mm)	10 pcs
REF P6850	ParaPost XT Titanium Posts (.050"/1.25mm)	10 pcs
REF P6855	ParaPost XT Titanium Posts (.055"/1.40mm)	10 pcs
REF P6860	ParaPost XT Titanium Posts (.060"/1.50mm)	10 pcs
REF P7433	ParaPost XP Plastic Impression Posts (.036"/0.90mm)	20 pcs
REF P7434	ParaPost XP Plastic Impression Posts (.040"/1.00mm)	20 pcs
REF P74345	ParaPost XP Plastic Impression Posts (.045"/1.14mm)	20 pcs
REF P7435	ParaPost XP Plastic Impression Posts (.050"/1.25 mm)	20 pcs
REF P74355	ParaPost XP Plastic Impression Posts (.055"/1.40 mm)	20 pcs
REF P7436	ParaPost XP Plastic Impression Posts (.060"/1.50 mm)	20 pcs
REF P7437	ParaPost XP Plastic Impression Posts (.070"/1.75mm)	20 pcs
REF P7463	ParaPost XP Titanium Temporary Posts (.036"/0.90mm)	20 pcs
REF P7464	ParaPost XP Titanium Temporary Posts (.040"/1.00mm)	20 pcs
REF P74645	ParaPost XP Titanium Temporary Posts (.045"/1.14mm)	20 pcs
REF P7465	ParaPost XP Titanium Temporary Posts (.050"/1.25mm)	20 pcs
REF P74655	ParaPost XP Titanium Temporary Posts (.055"/1.40mm)	20 pcs
REF P7466	ParaPost XP Titanium Temporary Posts (.060"/1.50mm)	20 pcs
REF P7467	ParaPost XP Titanium Temporary Posts (.070"/1.75mm)	20 pcs
REF P7513	ParaPost XP Plastic Burnout Posts (.036"/0.90mm)	10 pcs
REF P7514	ParaPost XP Plastic Burnout Posts (.040"/1.00mm)	10 pcs
REF P75145	ParaPost XP Plastic Burnout Posts (.045"/1.14mm)	10 pcs
REF P7515	ParaPost XP Plastic Burnout Posts (.050"/1.25mm)	10 pcs
REF P75155	ParaPost XP Plastic Burnout Posts (.055"/1.40mm)	10 pcs
REF P7516	ParaPost XP Plastic Burnout Posts (.060"/1.50mm)	10 pcs
REF P7517	ParaPost XP Plastic Burnout Posts (.070"/1.75mm)	10 pcs
REF P7843	ParaPost XP Titanium Posts (.036"/0.90mm)	10 pcs
REF P7844	ParaPost XP Titanium Posts (.040"/1.00mm)	10 pcs
REF P78445	ParaPost XP Titanium Posts (.045"/1.14mm)	10 pcs
REF P7845	ParaPost XP Titanium Posts (.050"/1.25mm)	10 pcs
REF P78455	ParaPost XP Titanium Posts (.055"/1.40mm)	10 pcs
REF P7846	ParaPost XP Titanium Posts (.060"/1.50mm)	10 pcs
REF P884	ParaPost XH Titanium Posts (.040"/1.00mm)	10 pcs
REF P8845	ParaPost XH Titanium Posts (.045"/1.14mm)	10 pcs
REF P883	ParaPost XH Titanium Posts (.036"/0.90mm)	10 pcs
REF P885	ParaPost XH Titanium Posts (.050"/1.25mm)	10 pcs
REF P8855	ParaPost XH Titanium Posts (.055"/1.40mm)	10 pcs
REF P886	ParaPost XH Titanium Posts (.060"/1.50mm)	10 pcs
REF P887	ParaPost XH Titanium Posts (.070"/1.75mm)	10 pcs



OFFER

**3+1\***

ON REFILLS





- Global market leader in post systems
- Proven clinical success with more than 50 studies
- More than 55 years of expertise
- One-office-visit and laboratory techniques
- Endo meets Resto – complete system with core-build ups and cements

## ParaPost System Kits

REF P780T	XP Intro Kit One-Office-Visit Technique Titanium Alloy
REF P880	XH Intro Kit One-Office-Visit Titanium
REF P780	XP Intro Kit One-Office-Visit Stainless Steel
REF P781	XP Casting Technique Intro Kit
REF P680T	XT One-Office-Visit Intro Kit Titanium
REF PF170	Fiber Lux Intro Kit
REF PF160	Fiber White Intro Kit
REF PF180	TaperLux Intro Kit/One-Office-Visit Technique



**Fiber Lux Introductory Kit**  
Metal-free Posts  
REF PF170



**X System Introductory Kit**  
REF P780T XP Titanium posts  
REF P780 XP Stainless steel posts  
REF P880 XH Titanium posts  
REF P680T XT Titanium posts



**Taper Lux Introductory Kit**  
Metal-free Posts  
REF PF180



**System Introductory Kit**  
Casting  
REF P781 XP



LEARN MORE



# A-silicone impression material

## AFFINIS



Still one of the fastest A-silicones



Precise moulding of every detail even under the most difficult conditions



Confidence during impression taking



Easy monitoring for success



### AFFINIS putty

2 x 300ml Base and Catalyst  
2 x Dosage spoon

#### Viscosity

putty soft		6530
putty soft fast		6531
putty super soft		6535



### AFFINIS system 360

2 x 380ml Cartridge  
1 x Fixation Ring

#### Viscosity





putty		6472
heavy body		6487
MonoBody		6482



### AFFINIS and AFFINIS PRECIOUS

2 x 50ml Cartridge  
12 x Universal Mixing Tip

#### Viscosity

light		6501
light fast		6601
regular		6511
regular fast		6611
PRECIOUS light		6775
PRECIOUS regular		6776

**OFFER**

**3+1\***

ON REFILLS

**FOR THE FULL RANGE**  
PLEASE ASK YOUR TERRITORY MANAGER

LEARN MORE





A-SILICONE IMPRESSION MATERIAL

# PRESIDENT

The Original

- High tear strength and elasticity
- Optimised flowing behavior and thixotropy
- Easier-to-mix putties
- New and more contrasting colours

**OFFER**



**3+1\***

ON REFILLS



## PRESIDENT The Original putty

2 x 300ml Base and Catalyst  
2 x Spoons

Viscosity		
putty		60019942
putty soft		60019943
putty soft fast		60022166
putty super soft		60019944



## PRESIDENT The Original system 360

2 x 380ml Cartridge  
1 x Fixation Ring


Viscosity		
MonoBody		60019946
heavy body		60019948

**FOR THE FULL RANGE**  
PLEASE ASK YOUR TERRITORY MANAGER



## PRESIDENT The Original light cartridge 50ml

2 x 50ml Cartridge  
12 x Universal Mixing Tips

Viscosity		
regular body		60019939
light body		60019938
Xtra light body		60019937

LEARN MORE





Not just any wipe!

# SciCan OPTIM 1



ONE-STEP CLEANER  
& DISINFECTANT

Active Ingredient: 0.5% Hydrogen Peroxide



- Disinfection in just one minute
- Smear free cleaning
- Durable and wet for covering a wide surface area
- Safe to use on upholstery and hard surfaces
- Recommended for intra oral cameras and scanners



**SciCan OPTIM 1  
Ready-to-use  
Pre-saturated Wipes**

**Standard Case**  
Code: OPTEU1-W12  
12 standard tubs of 160 wipes

**Standard Tub,**  
160 wipes, 15cm x 18cm

ASK YOUR  
TERRITORY  
MANAGER FOR A  
DEMONSTRATION

LEARN MORE



\*\* Effective within 1 minute according to European Standards; EN 14476, EN 14348, EN 13727, EN 13624, EN 16615. According to the required RKI/DVV standard in Germany SciCan OPTIM 1 is FULL VIRUCIDAL within 2 minutes





## Ultrasonic Cleaning System

Ultrasonic cleaning supports safe instrument reprocessing!

# BioSonic UC150

- 5.7l stainless steel tank
- Low noise during operation
- Data log function and download of usage data via USB
- Individually adjustable functions
- Offers space for several cassettes at the same time
- Use as tabletop or built-in unit
- Heating function from 25 up to 60 °C (Function can be activated separately via code)
- Glove-friendly touch screen



LEARN MORE



### BioSonic UC150 Cleaner Unit UK Plug type G

REF 60021989

- 1 × Unit with display
- 1 × Recess Template
- 1 × Basket
- 1 × Recess Gasket
- 1 × Power Cord [UK]
- 1 × Drain Tube

### Accessories

- REF UC53 1 × 600 ml Beaker with Cover & Positioning Ring
- REF 60022040 1 × Cassette Basket 26.8 × 15.8 × 11.5 cm
- REF 60022041 1 × Beaker Holder

### Ultrasonic Solutions

- REF 60021960 1 × BioSonic UC40 dosed bottle, 1l
- REF 60021961 1 × BioSonic UC40 refill can, 5l



# BioSonic UC50DB

### BioSonic UC50DB Cleaner Unit UK Plug

REF UC50DB230UK

- With Digital Timer
- Tank capacity 3.2 l
- 1 × Tank cover
- 1 × Fingerguard instrument basket
- 1 × Owners guide



LEARN MORE



NOW IN NEW BLISTER PACKS

DIATECH

# Diamond Burs

DIATECH Diamonds now come in individually sealed and marked blister packs for traceability and convenience.

- Gold-plated premium instruments
- Higher number of diamonds per mm<sup>2</sup>
- 100% natural diamond grit
- 100% visual quality control
- ISO standardized grit sizes
- Precise blank calibration

**\*Large Size**

Content: 50 instruments  
 10 × 5 blister packed burs  
 Available for selected products

**Regular Size**

Content: 5 instruments  
 Packed in 5 single sealed blister trays  
 Available for all products



NOW AVAILABLE IN PACKS OF **50\***

LEARN MORE



**Order Information**

Head shape

Head shape number

Head diameter in 1/10 mm

Head length in mm

Grit sizes:

C	107–181 μm*
M	64–126 μm*
F	27–76 μm*
EF	10–36 μm*

\* Grit sizes according to ISO 7711-3

Head shape				Round End Taper			
Price Gr.	C	M	F	EF	D1/005D	D1/005D	D1/005D
Head shape no.					856		
ISO Ø 1/10 mm					18		
Head L in mm					9		
					60032312	60032318	60032316
							Y
							Z
							60032314

Illustration

Y: This item is also available with friction grip short shank (FGSS)

Z: This item is also available in packaging units of 50 pieces.

For article numbers see table below.


































































































Price Group



# DIATECH Diamond Bur Range

Head shape	Head shape no.	ISO Ø 1/10 mm	Head L in mm	Price Gr.	C	M	F	EF
Bud/Pointed Football	368	16	3.5	D1/00SD	60031758			
Bud/Pointed Football	368	20	5	D1/00SD	60031756			
Bud/Pointed Football	368	23	5	D1/00SD	60031754			
Bud/Pointed Football	368	23	5	D6/00SD	60031760			
Bud/Pointed Football	368	23	5	D1/00SD	60031766			
Bud/Pointed Football	368	23	5	D1/00SD	60031762			
Bud/Pointed Football	368	23	5	D6/00SD	60031770			
Bud/Pointed Football	368	23	5	D1/00SD	60031776			
Bud/Pointed Football	368	23	5	D1/00SD	60031772			
Egg/Football	379	16	3.5	D1/00SD	60031786			
Egg/Football	379	18	3.5	D1/00SD	60031784			
Egg/Football	379	18	3.5	D1/00SD	60031782			
Egg/Football	379	23	5	D1/00SD	60031794			
Egg/Football	379	23	5	D6/00SD	60031792			
Egg/Football	379	23	5	D6/00SD	60031788			
Egg/Football	379	23	5	D6/00SD	60031796			
Egg/Football	379	23	5	D1/00SD	60031810			
Egg/Football	379	23	5	D1/00SD	60031805			
Egg/Football	379	23	5	D1/00SD	60031800			
Bullet	390	14	3	D1/00SD	60031818			
Bullet	390	14	3	D1/00SD	60031816			
Interproximal	392	16	6	D6/00SD	60031826			
Interproximal	392	16	6	D6/00SD	60031824			
Interproximal	392	16	6	D6/00SD	60031820			
Interproximal	392	18	8	D1/00SD	60031830			
Interproximal	392	18	8	D1/00SD	60031828			
Interproximal	392	18	8	D6/00SD	60031836			
Interproximal	392	18	8	D6/00SD	60031834			
Interproximal	392	21	10	D6/00SD	60031832			
Interproximal	392	21	10	D6/00SD	60031840			
Interproximal	392	21	10	D6/00SD	60031838			
Round	801	9	0.9	D1/00SD	60031850			
Round	801	21	1	D1/00SD	60031852			
Diabolic/Amalgam Remover	813	18	2	D6/00SD	60031952			
Wheel	818	25	0.7	D5/00SD	60031954			
Wheel	818	50	0.7	D5/00SD	60031956			
Wheel	828	22	1	D5/00SD	60031958			
Wheel	828	26	1	D5/00SD	60031960			
Wheel	828	30	1	D5/00SD	60031962			
Pear	830	9	2.7	D6/00SD	60031966			
Pear	830	10	2.7	D1/00SD	60031968			
Pear	830	12	2.7	D1/00SD	60031970			
Pear	830	14	2.7	D1/00SD	60031972			
Pear	830	16	2.7	D1/00SD	60031974			
Depth Marker	834	16	6.8	D5/00SD	60031996			
Depth Marker	834	21	6.8	D5/00SD	60031998			
Flat End Taper	847	14	8	D1/00SD	60032130			
Flat End Taper	847	16	8	D6/00SD	60032128			
Flat End Taper	847	16	8	D6/00SD	60032132			
Flat End Taper	847	18	8	D1/00SD	60032136			
Flat End Taper	847	18	8	D1/00SD	60032134			
Flat End Taper	847	18	8	D6/00SD	60032138			
Flat End Taper	847	25	8	D6/00SD	60032142			
Flat End Taper	847	25	8	D6/00SD	60032140			
Flat End Taper	847	25	8	D6/00SD	60032144			
Flat End Taper	847	25	8	D6/00SD	60032146			
Flat End Taper	848	12	10	D1/00SD	60032160			
Flat End Taper	848	14	10	D1/00SD	60032164			
Flat End Taper	848	14	10	D1/00SD	60032162			
Flat End Taper	848	16	10	D6/00SD	60032166			
Flat End Taper	848	16	10	D1/00SD	60032170			
Flat End Taper	848	16	10	D1/00SD	60032168			
Flat End Taper	848	18	10	D6/00SD	60032172			
Flat End Taper	848	18	10	D1/00SD	60032174			
Flat End Taper	848	23	10	D6/00SD	60032176			
Flat End Taper	848	23	10	D6/00SD	60032180			
Flat End Taper	848	23	10	D6/00SD	60032178			
Round End Taper	849	16	4	D1/00SD	60032190			
Round End Taper	849	25	4	D6/00SD	60032192			
Round End Taper	849	25	4	D6/00SD	60032194			
Round End Taper	850	12	10	D1/00SD	60032196			
Round End Taper	850	14	10	D1/00SD	60032200			
Round End Taper	850	14	10	D1/00SD	60032204			
Round End Taper	850	14	10	D1/00SD	60032202			
Needle	859	12	10	D1/00SD	60032408			
Needle	859	14	10	D1/00SD	60032416			
Needle	859	14	10	D1/00SD	60032413			
Needle	859	16	10	D1/00SD	60032410			
Needle	859	16	10	D1/00SD	60032418			
Needle	859	18	10	D1/00SD	60032424			
Needle	859	18	10	D1/00SD	60032422			
Needle	859	18	10	D6/00SD	60032420			
Flame	860	10	4	D6/00SD	60032426			
Flame	860	12	5	D6/00SD	60032432			
Flame	860	12	5	D6/00SD	60032430			
Flame	860	12	5	D6/00SD	60032428			
Flame	860	14	5	D6/00SD	60032442			
Flame	860	14	5	D6/00SD	60032440			
Flame	860	14	5	D6/00SD	60032438			
Flame	862	12	8	D1/00SD	60032450			
Flame	862	12	8	D1/00SD	60032448			
Flame	862	12	8	D1/00SD	60032446			
Flame	862	14	8	D1/00SD	60032456			
Flame	862	14	8	D1/00SD	60032454			
Flame	862	14	8	D1/00SD	60032452			
Flame	862	16	8	D1/00SD	60032460			
Flame	862	16	8	D1/00SD	60032458			
Flame	862	16	8	D1/00SD	60032466			
Flame	862	16	8	D1/00SD	60032464			
Flame	862	16	8	D1/00SD	60032462			
Flame	863	12	10	D1/00SD	60032474			
Flame	863	12	10	D1/00SD	60032472			
Flame	863	12	10	D1/00SD	60032470			
Acorn	905	27	2.9	D4/00SD	60032668			
Acorn	905	31	3.1	D4/00SD	60032670			
Round Wheel	909	40	1	D5/00SD	60032672			
Round Wheel	909	40	1.5	D5/00SD	60032674			
Round Wheel	909	40	1.5	D5/00SD	60032676			
Round Wheel	909	55	2	D5/00SD	60032678			
Strated Cylinder	6051	18	6	D2/00SD	60031842			
Strated Cylinder	6052	18	8	D2/00SD	60031844			
Strated Cylinder	6052	18	8	D2/00SD	60031846			
Strated Cylinder	6055	18	8	D2/00SD	60031848			
Pear Long	830L	10	4	D1/00SD	60031976			
Pear Long	830L	12	4	D1/00SD	60031982			
Pear Long	830L	14	5	D6/00SD	60031980			
Pear Long	830L	14	5	D6/00SD	60031978			
Pear Long	830L	16	5	D1/00SD	60031986			
Pear Long	830L	16	5	D6/00SD	60031984			
Pear Long	830L	16	5	D1/00SD	60031988			
Pear Long	830L	16	5	D1/00SD	60031990			
Pear Long	830L	18	5	D1/00SD	60031992			
Round End Taper with Pin	856P	18	8	D3/00SD	60032364			
Round End Taper with Pin	856P	21	8	D3/00SD	60032366			
Round End Taper Extra Long	856XL	16	12	D1/00SD	60032376			
Round End Taper Extra Long	856XL	18	12	D1/00SD	60032378			
Tapered Torpedo /Curette	877K	12	6	D6/00SD	60032493			
Tapered Torpedo /Curette	877K	14	6	D6/00SD	60032496			
Tapered Torpedo /Curette	877K	16	6	D6/00SD	60032498			
Tapered Torpedo /Curette	877K	18	6	D6/00SD	60032500			
Tapered Torpedo /Curette	878K	10	8	D6/00SD	60032516			
Tapered Torpedo /Curette	878K	12	8	D6/00SD	60032518			
Tapered Torpedo /Curette	878K	14	8	D6/00SD	60032524			
Tapered Torpedo /Curette	878K	14	8	D6/00SD	60032522			
Tapered Torpedo /Curette	878K	16	8	D6/00SD	60032526			
Tapered Torpedo /Curette	878K	16	8	D6/00SD	60032530			
Tapered Torpedo /Curette	878K	18	8	D6/00SD	60032528			
Tapered Torpedo /Curette	878K	18	8	D6/00SD	60032532			
Tapered Torpedo /Curette	878K	18	8	D6/00SD	60032536			
Tapered Torpedo /Curette	878K	18	8	D6/00SD	60032534			
Flat End Cylinder	837T	14	8	D3/00SD	60032070			
Flat End Cylinder	837T	16	8	D3/00SD	60032072			
Flat End Taper	847T	16	8	D3/00SD	60032136			
Flat End Taper	847T	18	8	D3/00SD	60032138			
Flat End Taper	848T	16	10	D3/00SD	60032186			
Flat End Taper	848T	18	10	D3/00SD	60032188			
Round End Taper	850T	14	10	D3/00SD	60032234			
Round End Taper	850T	16	10	D3/00SD	60032236			
Round End Taper	850T	18	10	D3/00SD	60032238			
Round End Taper	850T	18	10	D3/				

# DIATECH Diamond Bur Range (cont.)

														
	801 12 1.2	801 14 1.4	801 16 1.6	801 18 1.8	801 20 2	801 23 2.3	801 27 2.7	801 29 2.9	802 12 1.2	802 14 1.4	802 16 1.6	805 12 0.9	805 14 1.4	805 16 1.6
	D1/005D 60031856	D1/005D 60031862	D6/005D 60031858	D1/005D 60031867	D6/005D 60031864	D1/005D 60031875	D1/005D 60031872	D6/005D 60031870	D1/005D 60031880	D6/005D 60031878	D6/005D 60031882	D1/005D 60031889	D1/005D 60031886	D6/005D 60031884
														
	835 8 3	835 10 4	835 12 4	835 14 4	835 16 4	835 18 4	836 12 6	836 14 6	837 12 7	837 12 8	837 14 8	837 16 8	838 9 3	838 10 4
	D1/005D 60032002	D1/005D 60032005	D1/005D 60032008	D1/005D 60032010	D1/005D 60032012	D1/005D 60032028	D6/005D 60032026	D6/005D 60032030	D1/005D 60032032	D1/005D 60032050	D1/005D 60032052	D1/005D 60032054	D6/005D 60032056	D1/005D 60032058
														
	850 15 10	850 15 10	850 16 10	850 18 10	850 23 10	851 8 8	851 8 8	852 12 6	852 16 6	852 16 6	853 8 3.5	855 10 6	855 12 6	856 14 8
	D6/005D 60032206	D1/005D 60032210	D1/005D 60032208	D6/005D 60032212	D1/005D 60032215	D6/005D 60032218	D1/005D 60032222	D1/005D 60032220	D6/005D 60032224	D1/005D 60032246	D1/005D 60032248	D1/005D 60032248	D1/005D 60032270	D1/005D 60032272
														
	863 16 10	864 14 12	877 10 6	877 12 6	878 10 8	878 12 8	878 14 8	878 16 8	879 12 10	879 14 10	879 16 10	879 16 10	880 12 6	881 14 6
	D1/005D 60032480	D1/005D 60032478	D1/005D 60032476	D1/005D 60032486	D1/005D 60032484	D6/005D 60032488	D6/005D 60032490	D1/005D 60032504	D1/005D 60032502	D1/005D 60032508	D1/005D 60032506	D1/005D 60032510	D1/005D 60032512	D1/005D 60032554
														
	830L 20 5	835R 8 3	835R 10 4	835R 12 4	835R 14 4	836R 8 6	836R 10 6	836R 12 6	836R 14 6	836R 14 6	837R 12 8	837R 14 8	837R 16 8	845R 16 4
	D1/005D 60031994	D6/005D 60032016	D1/005D 60032018	D1/005D 60032022	D6/005D 60032020	D6/005D 60032024	D1/005D 60032034	D1/005D 60032038	D1/005D 60032036	D1/005D 60032040	D1/005D 60032044	D1/005D 60032042	D1/005D 60032062	D1/005D 60032060
														
	878K 21 8	878K 18 8	878K 21 8	879K 12 10	879K 14 10	879K 16 10	879K 18 10	879K 21 10	879K 23 10	879K 23 10	801M 9 0.9	838M 7 2.5	889M 7 2.5	ET3 8 3
	D6/005D 60032538	D6/005D 60032540	D3/005D 60032544	D3/005D 60032542	D3/005D 60032548	D3/005D 60032546	D6/005D 60032564	D6/005D 60032566	D6/005D 60032572	D6/005D 60032570	D6/005D 60032576	D6/005D 60032574	D6/005D 60032578	D6/005D 60032580
														
	856T 18 8	856T 21 9	856T 25 7	878KT 16 8	879KT 14 10	879KT 16 10	881T 14 8	881T 16 8	882T 14 10	882T 14 10	379Z 23 5	801Z 12 1.2	801Z 14 1.4	801Z 18 1.8
	D3/005D 60032370	D3/005D 60032372	D3/005D 60032374	D3/005D 60032550	D3/005D 60032584	D3/005D 60032586	D3/005D 60032616	D3/005D 60032618	D3/005D 60032630	D5/005D 60031814	D5/005D 60031900	D5/005D 60031902	D5/005D 60031904	D5/005D 60032242





Inverted Cone with Collar	806 10 2.2	D6/005D 60031918	Inverted Cone with Collar	806 12 2.2	D6/005D 60031920	Inverted Cone with Collar	806 14 2.2	D6/005D 60031922	Inverted Cone with Collar	806 16 2.5	D6/005D 60031924	Inverted Cone	807 10 3.5	D6/005D 60031926	Inverted Cone	807 12 3.5	D6/005D 60031928	Inverted Cone	807 16 4	D6/005D 60031930	Inverted Cone	807 18 5	D6/005D 60031932	Occlusal Reduction/ Barrel	811 31 4.2	D4/005D 60031934	Occlusal Reduction/ Barrel	811 33 5	D4/005D 60031936	Occlusal Reduction/ Barrel	811 37 7	D4/005D 60031940	Occlusal Reduction/ Barrel	811 48 8	D4/005D 60031942	Occlusal Reduction/ Barrel	811 58 10	D4/005D 60031944	Diabolo/Amalgam Remover	813 15 1.5	D6/005D 60031948	Diabolo/Amalgam Remover	813 14 1.5	D6/005D 60031950																																										
Round End Cylinder	838 14 4	D1/005D 60032080	End Cutting Cylinder	839 10 0.2	D2/005D 60032084	End Cutting Cylinder	839 12 0.2	D2/005D 60032088	End Cutting Cylinder	839 14 0.2	D2/005D 60032086	End Cutting Cylinder	839 16 0.2	D2/005D 60032092	End Cutting Cylinder	839 18 0.2	D2/005D 60032090	Flat End Cylinder	842 14 12	D1/005D 60032094	Flat End Taper	845 10 4	D6/005D 60032096	Flat End Taper	845 12 4	D6/005D 60032098	Flat End Taper	845 14 4	D6/005D 60032100	Flat End Taper	845 16 4	D6/005D 60032104	Flat End Taper	845 18 4	D6/005D 60032102	Flat End Taper	846 16 6	D1/005D 60032112	Flat End Taper	846 12 7	D1/005D 60032114	Flat End Taper	846 18 7	D1/005D 60032116	Flat End Taper	846 23 7	D6/005D 60032118	Flat End Taper	847 12 8	D1/005D 60032126																																				
Round End Taper	856 14 9	D1/005D 60032286	Round End Taper	856 16 9	D6/005D 60032288	Round End Taper	856 18 9	D1/005D 60032292	Round End Taper	856 20 9	D1/005D 60032290	Round End Taper	856 22 9	D6/005D 60032300	Round End Taper	856 24 9	D1/005D 60032304	Round End Taper	856 26 9	D6/005D 60032302	Round End Taper	856 28 9	D6/005D 60032306	Round End Taper	856 30 9	D1/005D 60032310	Round End Taper	856 32 9	D1/005D 60032314	Round End Taper	856 34 9	D6/005D 60032318	Round End Taper	856 36 9	D1/005D 60032316	Round End Taper	856 38 9	D6/005D 60032320	Round End Taper	856 40 9	D6/005D 60032322	Round End Taper	856 42 9	D6/005D 60032326	Round End Taper	856 44 9	D6/005D 60032324	Round End Taper	856 46 9	D6/005D 60032328	Round End Taper	856 48 9	D6/005D 60032332	Round End Taper	856 50 9	D6/005D 60032330	Needle	858 14 8	D1/005D 60032388	Needle	858 16 8	D1/005D 60032386	Needle	858 18 8	D1/005D 60032384	Needle	858 20 8	D1/005D 60032394	Needle	858 22 8	D1/005D 60032392	Needle	858 24 8	D1/005D 60032390	Needle	858 26 8	D1/005D 60032396	Needle	858 28 8	D1/005D 60032400	Needle	858 30 8	D1/005D 60032401	Needle	858 32 8	D1/005D 60032398
Round End Cylinder	881 12 8	D6/005D 60032598	Round End Cylinder	881 14 8	D1/005D 60032602	Round End Cylinder	881 16 8	D1/005D 60032600	Round End Cylinder	881 18 8	D6/005D 60032604	Round End Cylinder	881 20 8	D1/005D 60032608	Round End Cylinder	881 22 8	D1/005D 60032606	Round End Cylinder	881 24 8	D6/005D 60032610	Round End Cylinder	881 26 8	D1/005D 60032614	Round End Cylinder	881 28 8	D1/005D 60032612	Beveled Cylinder	882 12 6	D6/005D 60032622	Beveled Cylinder	882 14 6	D6/005D 60032620	Beveled Cylinder	882 16 6	D6/005D 60032624	Beveled Cylinder	882 18 6	D6/005D 60032628	Beveled Cylinder	882 20 6	D6/005D 60032632	Beveled Cylinder	882 22 6	D6/005D 60032636	Beveled Cylinder	882 24 6	D6/005D 60032640	Beveled Cylinder	882 26 6	D6/005D 60032642	Beveled Cylinder	882 28 6	D6/005D 60032646	Beveled Cylinder	882 30 6	D6/005D 60032644	Beveled Cylinder	882 32 6	D6/005D 60032648	Beveled Cylinder	882 34 6	D6/005D 60032652	Beveled Cylinder	882 36 6	D6/005D 60032654	Beveled Cylinder	882 38 6	D6/005D 60032658	Beveled Cylinder	882 40 6	D6/005D 60032656	Flange	889 9 3.5	D1/005D 60032660	Occlusal Reduction	889 12 7	D4/005D 60032664	Acorn	905 23 2.7	D4/005D 60032666						
Rounded Edge Taper	846R 16 6	D1/005D 60032122	Rounded Edge Taper	847R 16 8	D6/005D 60032148	Rounded Edge Taper	847R 18 8	D1/005D 60032152	Rounded Edge Taper	847R 20 8	D1/005D 60032150	Rounded Edge Taper	847R 22 8	D6/005D 60032154	Rounded Edge Taper	848R 16 10	D6/005D 60032182	Rounded Edge Taper	848R 18 10	D6/005D 60032184	Round End Taper Long	850L 16 12	D1/005D 60032228	Round End Taper Long	850L 18 12	D1/005D 60032230	Round End Taper Narrow	850N 14 10	D6/005D 60032232	Round End Taper Long	856L 12 10	D1/005D 60032340	Round End Taper Long	856L 14 10	D1/005D 60032342	Round End Taper Long	856L 16 10	D6/005D 60032344	Round End Taper Long	856L 18 10	D1/005D 60032350	Round End Taper Long	856L 20 11	D1/005D 60032348	Round End Taper Long	856L 22 11	D6/005D 60032352	Round End Taper Long	856L 24 11	D1/005D 60032356	Round End Taper Long	856L 26 11	D1/005D 60032354	Round End Taper Narrow	856N 14 8	D6/005D 60032362																														
ET4	10 4.2	D2/005D 60032713	ET6	14 6	D2/005D 60032715	ET9	14 9	D2/005D 60032714	KS0	10 7	D2/005D 60032698	KS1	12 7	D2/005D 60032717	KS2	14 8	D2/005D 60032718	KS4	24 4.2	D2/005D 60032719	KS5	31 6.2	D2/005D 60032721	KS6	25 9	D2/005D 60032722	KS7	54 1	D2/005D 60032723	Egg/Football	379T 5	D3/005D 60031812	Occlusal Reduction/Barrel	811T 33 5	D3/005D 60031946	Flat End Cylinder	836T 14 6	D3/005D 60032046	Flat End Cylinder	836T 16 6	D3/005D 60032048																																													
Round End Taper	850Z 18 10	D5/005D 60032444	Round End Taper	856Z 18 8	D5/005D 60032380	Flame	862Z 14 8	D5/005D 60032468	Flame	863Z 14 10	D5/005D 60032482	Round End Cylinder	881Z 14 8	D5/005D 60032620	Round End Cylinder	881Z 16 8	D5/005D 60032623																																																																					

# MEET THE TEAM



**Laura Whyte**  
07788 146 109  
laura.whyte@coltene.com  
Scotland



**Debi Boynton**  
07764 327 425  
debi.boynton@coltene.com  
North West



**Sue Robins**  
07786 965 019  
sue.robins@coltene.com  
East Midlands




**Emily Downie**  
07889 227 424  
emily.downie@coltene.com  
North East



**Nick O'Keeffe**  
00353 8641 35766  
nicholas.okeeffe@coltene.com  
Ireland



**Emma Gilks**  
07906 749 682  
emma.gilks@coltene.com  
West Midlands & North Wales




**Shah Monnan**  
07712 407 071  
shah.monnan@coltene.com  
Essex & East London



**Andrew Miller**  
07788 583 057  
andrew.miller@coltene.com  
South Wales & West Country



**Eva Fedorowicz**  
07818 424 210  
eva.fedorowicz@coltene.com  
North, West & Central London



**Louise Fish**  
07788 146 130  
louise.fish@coltene.com  
South London & South Coast

## UK OFFICE - SUSSEX



**Nic Coomber**  
01444 235 486 - Ext 204  
National Account Manager  
& Marketing Strategist



**Helen Chart**  
01444 235 486 - Ext 202  
Finance & Office Manager

\*Offer ends 31.12.2024. Available from participating suppliers. RRP's may vary with supplier.  
For free item offers, free item from supplier unless otherwise stated.

© COLTENE – Printing errors, mistakes and changes reserved.



0800 254 5115 | [info.uk@coltene.com](mailto:info.uk@coltene.com) | [www.coltene.com](http://www.coltene.com)

Coltene/Whaledent Ltd., The President Suite A, Kendal House, Victoria Way, Burgess Hill, West Sussex RH15 9NF

