



Oro Clean Chemie AG  
Allmendstrasse 21  
8320 Fehraltorf

**Efficacy of DENTIRO® Sensitive against *Feline coronavirus*, strain NADL-2, RVB-1259 in a quantitative suspension test at 20 °C according to EN14476:2013+A1:2015 (E) under dirty condition**

**EXPERT OPINION**

This expert opinion is based on the test report VX-TR-18-0018 dated 10 January 2018.

The virucidal properties of the disinfectant DENTIRO® Sensitive of Oro Clean Chemie AG against *Feline coronavirus*, strain NADL-2, RVB-1259 was investigated by a quantitative suspension test according to EN14476:2013+A1:2015 (E) under dirty condition.

According to this suspension test, a disinfectant or a disinfectant solution at a particular concentration is considered as having virucidal activity if the virus titre is reduced by  $\geq 4 \log_{10}$  (inactivation  $\geq 99.99\%$ ) within the recommended exposure period.

DENTIRO® Sensitive was examined undiluted\* at 20 °C for the exposure time of 15 seconds and 30 seconds. After the exposure time, the virus reduction exceeded 4  $\log_{10}$ -steps in all assays. Therefore, a virucidal activity against *Feline coronavirus*, strain NADL-2, RVB-1259 was measured as follow:

Dirty condition	Undiluted*	15 seconds
Dirty condition	Undiluted*	30 seconds

Kuala Lumpur, 10 January 2018

**Suhana Ripim**  
Team Leader Virology Testing

**Zhao Min Khoo**  
Microbiologist

\* The product can only be tested at 80.00 % concentration or less, as some dilution always occurs when test organisms and interfering substance are added.