



COLTENE

April - June 2021

Promotions



CanalPro Jeni

→ p. 4 - 5



HyFlex Remover File

→ p. 7



DIATECH Instruments

→ p. 29



Sign up to the
COLTENE newsletter
and keep up to date.

CLICK OR SCAN

VacuSoft

NEW

Gentle and low-noise aspiration

Aspiration cannula distinguished by exceptional comfort, flexibility and silence!

- Soft tips provide gentle contact with the oral mucosa
- Low-noise aspiration without reflux
- Flexible adjustment to the treatment situation

VacuSoft tips are made of soft TPE plastic and are therefore particularly gentle in contact with the patient's oral mucosa, tongue and cheeks. The wide and flat cannula tip is ideally suited for keeping the patient's cheeks and tongue away and ensuring a clear view of the treatment field. Different degrees of hardness of the colour-coded tips provide flexibility in aspiration and allow for individual treatment situations.



VacuSoft sensitive comfort

Very soft tip for gentle aspiration in particularly (pressure-) sensitive patients



VacuSoft comfort

Soft tip for gentle aspiration in uncomplicated treatment cases



VacuSoft solid comfort

Stable tip for aspiration situations where more resistance is needed for retraction



3+1

BUY 3 REFILLS
GET 1 REFILL FREE*

Kit intro VacuSoft

REF 60023549

- 2 × VacuSoft Base
- 5 × Sensitive confort tip
- 5 × Comfort tip
- 5 × Solid tip



VacuSoft Base Refill

REF 60023550

- 10 × VacuSoft Base



VacuSoft Tip Refill

	50 ×	200 ×
Blue Tip Sensitive Comfort	REF 60023551	REF 60023556
Green Tip Comfort	REF 60023552	REF 60023557
Black Tip Solid Comfort	REF 60023554	REF 60023558



CanalPro™

Irrigation Solutions

Irrigation solutions are used to help dissolve vital and necrotic tissue, eradicate bacteria and other microorganisms, remove the smear layer, and lubricate the canal system.

As no single irrigation solutions exists that can accomplish all of the above goals, multiple solutions should be used during clinical treatment. Sodium hypochlorite (NaOCl), ethylenediaminetetraacetic acid (EDTA), and 2% Chlorhexidine (CHX) are commonly used for chemical debridement during the process of canal preparation.



CanalPro NaOCl

3% or 6%

For irrigation/debridement of root canals during and after treatment

Available in 3% and 6% formula

500ml - 3% - REF 60019655

500ml - 6% - REF 60019656



CanalPro EDTA 17%

17% EDTA solution (pH 8.5)

Removes smear layer and dentin mud

Opens dentin tubules for:

- Disinfection solutions to work better
- Better adhesion of sealers and obturation materials

500ml - REF 60019651



CanalPro CHX 2%

2% chlorhexidine digluconate solution

Substantive

Very effective against *E. faecalis* and fungal

500ml - REF 60019658

CanalPro™ Jeni

Digital assistance system for canal preparation

Autonomous "navigation" in the endo canal owing to revolutionary software

The new CanalPro Jeni endomotor has a patent-pending digital assistance system for canal preparation, which uses complex algorithms to control the file movement at millisecond intervals. The rotary motion and speed are continuously controlled based on current intensity and torque.

Different routes, obstacles and situations

CanalPro Jeni is named after its developer Prof. Dr. Eugenio Pedullà and guides you safely through the traffic like a navigation system. As such, the Jeni system automatically adapts to the individual root canal anatomy, working in unison with the mechanical and chemical preparation every step of the way.

Includes pre set programme settings for HyFlex and MicroMega rotary files

"The CanalPro Jeni is the best all-in-one endo motor and apex locator on the market.

As a dentist with practice limited to endodontics, I need reliable instruments, and the CanalPro Jeni is exactly that. This high quality machine has it all- accurate, crisp apex locator that responds to the file movements and gets very little interference from metal restorations, the motor responds immediately to wireless foot pedal, it has multiple file systems already installed and options to add more file systems, including reciprocation, if you wish. The motor adjusts itself according to the canal morphology, which helps reducing file separation.

I highly recommend this all-in-one endomotor for beginners and experienced alike."

Dr Sagi Shavit

DMD, MSc(Endodontology), MFGDP(UK), Dip.Dent.Imp



LEARN MORE



CLICK OR SCAN

Order a CanalPro Jeni motor and choose from one of the following free refill pack selections

4 x HyFlex CM

8 x HyFlex EDM

5 x MicroMega One Curve

Only available through your Territory Manager

Ask your Territory Manager for a demonstration and details of our new **LOYALTY SCHEME**

Jeni - ready - go!



CanalPro Jeni, complete set
REF 60023659

includes:

- 1 × Control unit
- 1 × Motor with apex measuring contact
- 1 × Fully insulated angled handpiece 1:1
- 1 × Bluetooth foot switch
- 1 × Power supply
- 1 × Apex cable set
- 1 × Holding bracket for apex cable

I use the files in combination with COLTENE's new CanalPro™ Jeni motor with integrated apex locator. The motor is smooth and efficient, and the audible signals provides the clinician with a regular reminder of the importance of constant irrigation. The Jeni has an integrated locator, useful as it provides information regarding the positioning of rotary files within the canals at all times, which could potentially reduce the incidence of overpreparation.

Dr Dhiraj Arora

BDS MJDF RCS(Eng) MSc (Endo) PG Cert CE



ASK ABOUT
OUR NEW
**JENI MOTOR
LOYALTY
SCHEME**

One Day Hands on Endodontic Course

with Dr Ammar Al Hourani

Dates

April	9	10
July	16	17
October	8	9
December	10	11

TO BOOK



CLICK OR SCAN

Venue Our Place, Unit 1, Windrose Point, Pynacles Close, Stanmore, HA7 4AG



Dr Ammar Al Hourani

BDS (Glasgow), DDSc (Endodontology), MFDS RCS Ed, MEndo RCS (Edin) Consultant/ Specialist in Endodontics and Clinical Lecturer, Peninsula Dental School

Ammar Al Hourani qualified from Glasgow Dental School in 2009 and gained an interest in Endodontics early on in his career. After several years as a Senior House Officer and General Dental Practitioner, his interest led him to embark on a 3-year full time postgraduate DDSc endodontic training programme at Liverpool Dental School in which he achieved his Membership in Endodontics specialism examination by the Royal College of Surgeons in Edinburgh (MEndo RCS Edin) and his Doctorate in Endodontology (DDSc).

Ammar is currently a Consultant and Clinical Lecturer at Peninsula Dental School, Plymouth, supervising the undergraduate endodontic teaching as well as, teaching on the Endodontic MSc programme at Peninsula. Alongside his university commitments, he currently teaches on several NHS 63 & DFT courses in the Yorkshire, Merseyside and East Midlands deaneries, as well as, eWisdom in London.

Ammar lectures nationally and internationally with many distinguished educational institutes and has a real passion for teaching and academia. This led him to create The EndoGuys, which provides endodontic teaching to GDPs across the UK.

A full day course where the delegates experience practical, safe and time saving root canal procedures using a combination of Micro Mega and Hyflex endodontic rotary files and to help move away from hand filing. Through hands on learning, the aim of the course is to give confidence and practical tips to overcome most of the endodontic cases where picking the right file for the procedure can reduce the risk of fracture. This in the long term will result in more confident practitioners, successful treatments and savings on surgery time and materials used in the practice. This is a practice builder for your business.

Delegates will work on extracted teeth/plastic blocks and use the latest rotary endodontic file systems from Coltene MicroMega and the Hyflex range.

- Learning objectives and CPD hours will be awarded as per GDC.
- Refreshments and lunch included.
- Socially distanced workstations with individual motors and materials plus PPE, ensure you have a safe and enjoyable learning experience.
- We supply the manufacturer's materials and equipment used on the course, plus a special loyalty scheme for all delegates who have attended the course.
- A personal referral network is created so your Dentists can discuss cases and predict possible outcomes on submitted cases.

This course provides your team with high quality training using quality materials from a dental manufacturer that believes your future clinical success is by upgrading your dentistry to a higher level.

Getting there Nearest tube station is Stanmore (Jubilee Line). Turn left out of the station - 10 mins walk. Nearest mainline station is Harrow & Wealdstone - approx 3 miles. Nearest free parking is Lady Aylesford Avenue - approx 8 mins walk. If you get lost on the day, call Kelly McCormack for help - **07818 424 210**

Day Agenda

- 8:30 - 9:00 Registration & Coffee
- 9:00 - 10:00 Modern Endodontic Workflow
- 10:15 - 10:30 Break
- 10:45 - 13:00 Hands on Training with Hyflex EDM and MicroMega OneCurve File Systems
- 13:00 - 14:00 Lunch
- 14:00 - 17:00 Hands on Training with Hyflex EDM and MicroMega OneCurve File Systems

GDC outcomes: C and D

Cost £350 inc. vat

Contact

info.uk@coltene.com
0800 254 5115

HyFlex™ Remover

NEW

The essential piece for your retreatment

One file to remove the canal obturation

Remover is a new addition to the Hyflex range for the removal of old obturation material, integrating the latest generation of cross-sectioning with 3 cutting edges and heat treatment.

Just one file for the removal of old material

You only need one instrument to remove the Gutta-Percha in the canal.

Efficient without solvent use


The only file on the market clinically qualified to be used without solvent.

Safe instrument

Due to its design, its non-cutting tip and its protocol, REMOVER avoids ledges and respects the root canal anatomy and while being effective.

Reduce the learning time

Remover fits perfectly into an existing retreatment procedure without requiring a change in approach. It's an essential tool for retreatment and can be used with conventional shaping files.



30/.07

Maximum torque : **2.5 N.cm**

Rotation speed : **400-800 rpm**

Heat treatment : **C.Wire**

HyFlex™ Remover

Blister pack of 5 instruments - Sterile

B5 HyFlex Remover L19 n30 7%

REF 60023649

B5 HyFlex Remover L23 n30 7%

REF 60023648



REMOVER FILE
IN ACTION



CLICK OR SCAN

"The Remover file has quickly turned into my preferred file for root canal retreatment. The tip of the file is sturdy enough to penetrate GP, and remove it quickly and effectively without the need of a solvent. Importantly, the file dimensions allow for a conservative preparation. I find it is flexible enough to work around canal curvatures and respects the canal anatomy well!"

Dr. Kreena Patel

BDS(Hons)MJDF
MClinDent(Dist) MEndo

Specialist Endodontist
Brigstock Dental Practice,
Thornton Heath, Surrey



ENDODONTICS





HyFlex™ EDM

STAYS ON TRACK

TREATS 80% OF ALL CASES WITH THE ONE FILE



HyFlex™ EDM Files

	Orifice Opener 25/.12	Glide Path File 10/.05	Preparation file 20/.05	OneFile 25/~	NEW FILE FOR COMPLEX CANALS Glide Path File 15/.03
					
	Refill × 3 Rotation speed 400 rpm and 2.5 Ncm	Refill × 3 Rotation speed 300 tr/min and 1.8 Ncm	Refill × 3 Rotation speed 400 rpm and 2.5 Ncm	Refill × 3 Rotation speed 400 rpm and 2.5 Ncm	Refill × 3
Long. 21 mm	-	REF 60022269	REF 60022275	REF 60022270	REF 60023493
Long. 25 mm	-	REF 60019594	REF 60022276	REF 60019595	REF 60023492
Long. 15 mm	REF 60019596	-	-	-	-

"I have been amazed at how predictably I can achieve excellent outcomes"

"Since my conversion to COLTENE HyFlex™ EDM files, I have been amazed at how predictably I can achieve excellent outcomes, even when taking on challenging teeth. I have gained confidence with the files' ability to effortlessly navigate sclerosed and curved canals. It is a simple system that can be used in almost all cases and it gives you a feeling of reassurance when you have such a high-quality file.

Their fracture resistance makes it possible for me to treat complex cases in my general practice setting, which is great for my patients.

Additionally, I have found the system to be very efficient, as it reduces the need for additional files, making endodontics more affordable".

Dr. Orry Torsney

Principal of Central Park Dental,
Mullingar, Co Westmeath.



6+2

BUY 6 FILE REFILLS
GET 2 REFILLS FREE*
MAY & JUNE ONLY

**ONLY AVAILABLE
THROUGH YOUR
COLTENE TERRITORY
MANAGER**

HyFlex™ EDM Glidepath Files

Glidepath HyFlex EDM File 15 / .03 GPF is a very secure solution for complex anatomies of the root canal.

- Easy instrumentation of very narrow canals
- Uncomplicated preparation with a limited field of vision in very narrow canals
- For creating a glide path in difficult anatomy



NEW FILE FOR COMPLEX CANALS

HyFlex EDM Glide Path File 15/.03, 21 mm
 REF 60023493
 3 x 15/.03

HyFlex EDM Glide Path File 15/.03, 25 mm
 REF 60023492
 3 x 15/.03

HyFlex™ EDM Finishing Files



HyFlex EDM Finishing File 40/.04
 REF 60019597

HyFlex EDM Finishing File 50/.03
 REF 60019598

HyFlex EDM Finishing File 60/.02
 REF 60019599

HyFlex™ EDM Kits

Simple shaping kit

REF 60022267 **21 mm**
 REF 60019363 **25 mm**

1 x Orifice Opener 25/12 

1 x Glide Path File 10/05 

1 x OneFile 25/~ 

NEW

Medium shaping kit

REF 60022277 **21 mm**
 REF 60022278 **25 mm**

1 x Glide Path File 10/05 

1 x Preparation File 20/05 

1 x OneFile 25/~ 

HyFlex EDM Finishing Kit (Optional)

REF 60019364

1 x 40/.04 

1 x 50/.03 

1 x 60/.02 

WATCH THE VIDEO




CLICK OR SCAN

HyFlex™ EDM Paper and Gutta-Percha Points

HyFlex paper and guttapercha points are optimally matched to the HyFlex EDM files and thus contribute to the best possible treatment result. They correspond to the respective sizes and have the perfect taper to dry and obturate root canals previously prepared with HyFlex files precisely and as quickly as possible.

HyFlex EDM guttapercha points are ideal for obturating root canals. HyFlex EDM paper points are ideal for drying root canals.

- Precise fit
- Colour-coded
- Highest quality
- assorted or individual sizes
- ≥ 28mm



HyFlex EDM Paper and Guttapercha points selection

Size / tape	Paper points 100 pcs	Guttapercha points 60 pcs
20 / .05	60022742	60022706
OneFile 25/~	60022743	60022707
40 / .04	60022676	60022708
50 / .03	60022741	60022709
60 / .02	60022740	60022710
Assortment 20/.05 - 60/.02	60022744	60022705



HyFlex™ CM GuttaPercha Points Greater Taper

HyFlex CM guttapercha points are ideal for obturating root canals which have been prepared with HyFlex CM files. At the tip of the point the diameter corresponds to the ISO sizes, but then the conicity follows the one from the CM files ≥ 28mm

HyFlex CM Guttapercha points Greater Taper selection

.04 #20	60 pcs	60022689
.06 #30	60 pcs	60022699
.06#45	60 pcs	60022702
.04 #60	60 pcs	60022696
.04 #20-45	60 pcs	60022703
.06 #20-45	100 pcs	60022704



Expansion meets adhesion

GuttaFlow[®] 2



GuttaFlow 2 is a filling system for obturation of root canals containing two products in one: Guttapercha in the form of powder with a particle size smaller than 30 µm and a sealer. This filling system works with cold, flowable Guttapercha.

- Excellent flow properties
- Almost no solubility
- Dense occlusion of the root canal
- Very good tissue compatibility
- Reinfection protection
- Optimal radiopacity
- Precise root post preparation possible

GuttaFlow[®] 2
Standard set 5ml syringe
 inc. 12 mixing tips
 REF 60013701

GuttaFlow[®] 2
Refill pack 5ml syringe
 REF 60013745

GuttaFlow[®] 2
Refill pack 24 mixing tips
 REF 60013719

Bioactive obturation material

GuttaFlow[®] bioseal




Not just for filling, also for healing!

- Reliable cold filling system for root canals
- Actively supports regeneration in the root canal
- Excellent flow properties
- Guttapercha and sealer in a single application
- Processing time: 5 minutes, curing time: 12 – 16 minutes

GuttaFlow[®] bioseal
Standard Set
 1 x Dual-barrel syringe,
 12 x Mixing tips,
 1 x Mixing pads
 REF 60019560

GuttaFlow[®] bioseal
Refill pack Applicator
 1 x Dual-barrel syringe
 REF 60019561

GuttaFlow[®] bioseal
Refill pack Mixing tips
 24 x Mixing tips
 REF 60019562



Images - Dr. Thomas Rieger

MicroMega One Curve

One single instrument for root canal shaping.
Heat treatment C.Wire

SMART

A single instrument for shaping

- ▶ Simplified management
- ▶ Single use: No risk of cross-contamination
- ▶ Reduced working time
- ▶ More time for Irrigation

EFFICIENT

C.Wire: Heat treatment of Nickel Titanium, A One Curve exclusive

- ▶ Easier access to the canal
- ▶ Pre-bendable
- ▶ Shape memory
- ▶ Conservation of the curvature



6+2

BUY 6 FILE REFILLS
- GET 2 REFILLS FREE*

MAY & JUNE ONLY

**ONLY AVAILABLE
THROUGH YOUR
COLTENE TERRITORY
MANAGER**



59% +
faster

Compared to an instrumental sequence ⁽¹⁾



33% +
faster

Compared to a heat-treated reciprocating instrument ⁽²⁾



2.4% X
more resistant
to fracture (cyclic fatigue)

100% +
more flexible



MicroMega One Curve


Blister pack of 5 instruments - Sterile

B5 One Curve n25 L25 4%	REF 20949441
B5 One Curve n35 L25 4%	REF 20949442
B5 One Curve n45 L25 4%	REF 20949443
B4 One Curve L25 Asst	REF 20949445
B5 One Curve n25 L31 4%	REF 20949451
B5 One Curve n25 L25 6%	REF 20949402
B5 One Curve n25 L25 6%	REF 20949412
Procedure pack	REF 20949412

One Curve Paper points and Guttapercha points

ONECURVE GP 25/.06 SB 60 pieces
REF 60023612

ONECURVE PP 25/.06 CP 144 pieces
REF 60023616



(1) Shaping ability of different single-file systems in severely curved root canals of extracted teeth - S. Burklein et al
(2) Data: MicroMega R&D

LEARN MORE



CLICK OR SCAN

REF 52002004

MicroMega Dual Move

- Cordless motor, well-balanced in hand and ambidextrous
- Miniature head, 340° rotative contra-angle
- Continuous rotation or reciprocation modes
- Fully customisable
- Accurate apex localisation
- Apex location – combined use

LEARN MORE



CLICK OR SCAN

SPECIAL OFFER PRICE *
SAVE 40%
ONLY £399 +vat

MicroMega Dual Pex*

- User friendly screen
- Accurate localisation even in wet and dry root canals, owing to multi-frequency positioning technology
- Stable delivered voltage whatever the battery level
- Real-time visual and sound indications
- Auto power-off to save the battery
- Plug & Play self tester

REF 51300050

* May vary with supplier

ParaPost® System

A complete range of posts for direct and indirect indications

- Global market leader in post systems
- Proven clinical success with more than 50 studies
- More than 55 years of expertise
- One-office-visit and laboratory techniques
- Endo meets Resto – complete system with core-build ups and cements



Buy a kit and receive a **SIZE 5 REFILL FREE***
(of the same kit type)

Fiber Lux Introductory Kit Metal-free Posts
REF PF170



Taper Lux Introductory Kit Metal-free Posts
REF PF180



System Introductory Kit Casting
REF P781 XP



X System Introductory Kit

- REF P780T XP Titanium posts
- REF P780 XP Stainless steel posts
- REF P880 XH Titanium posts
- REF P680T XT Titanium posts



LEARN MORE



CLICK OR SCAN

Dual Curing Composite Based Post Cement
and Core Build Up Material

ParaCore®

- 3 indications - 1 material
- For post cementation, core construction and cementation of crowns, bridges, inlays and onlays



ParaCore Refill Dentin REF 5887
2 × 5ml ParaCore Dentin
1 × 20 Mixing Tip Short Super Fine

ParaCore Refill White REF 5886
2 × 5ml ParaCore White
1 × 20 Mixing Tip Short Super Fine

ParaCore Slow Refill Dentin REF 60011391
2 × 5ml ParaCore Slow Dentin
1 × 20 Mixing Tip Short Super Fine

ParaCore Slow Refill Trans REF 60011392
2 × 5ml ParaCore Slow Trans
1 × 20 Mixing Tip Short Super Fine



ParaCore 5 ml Intro Kit REF 5885
1 × 5 ml ParaCore Dentin, 1 × 5 ml ParaCore White
1 × 5 ml ParaCore Dentin SLOW, 1 × 5 ml ParaCore Trans SLOW
3 × 3 ml ParaBond Non-Rinse, Conditioner + Adhesive A & B
20 × Mixing Tip Short Super Fine



LEARN MORE



CLICK OR SCAN



Flowable dual curing posterior restorative material

Fill-Up!®

Deep. Fast. Perfect.

- Optimal depth polymerisation with minimal shrinkage due to dual curing system.
- Guaranteed single-layer technique - even in very deep cavities.
- Optimised sealing of margins - reduced post-operative sensitivity
- Universal shade in a convenient Automix syringe for efficient placement



LEARN MORE



CLICK OR SCAN



Fill-Up! Intro Kit

REF 60019973

- 2 x 4.5 g Fill-Up!™ Universal
- 1 x 5 ml OneCoat 7 Universal
- 1 x 3 ml OneCoat 7 Universal Activator
- 6 x Mixing tips Short Super Fine
- 8 x Mixing tips Short Fine



Fill-Up! Refill

REF 60019341

- 1 x 4.5 g Fill-Up! Universal
- 3 x Mixing Tip Short Super Fine
- 4 x Mixing Tip Short Fine

" Given the move towards less amalgam use there is going to be greater need for suitable alternative materials and composite resin will lead the way. There are likely to be increased use of novel techniques including bulk-fill restorative materials. Fill-Up is one such material that allows rapid cavity filling, easy finishing and a reasonable wear resistance with simple top-up and repair options for long-term maintenance. "

Professor Brian Millar



The perfect way to try 3 new COLTENE materials

Restoration Highlight Kit

BRILLIANT EverGlow

- Exceptional polishability and long-lasting brilliance
- Aesthetic single-colour restorations
- Ideal handling through smooth consistency

ONE COAT 7 UNIVERSAL

- Self, Selective and Total Etch, a bond for all adhesive techniques
- Exceptional bonding values on enamel and dentine
- Excellent marginal seal and outstanding marginal leakage tightness

DIATECH ShapeGuard

- Adapts to any surface
- Uniform pressure distribution.
- Tooth morphology preserved
- Fast, easy and intuitive finishing & polishing
- Optimised for composite and ceramics



Restoration Highlight Kit Syringes

REF 60020090
 3 x BRILLIANT EverGlow Syringes (3g) (1 x A1/B1, A2/B2, A3/D3),
 1 x One Coat 7 Universal Bond 5ml +50 brushes,
 2 x DIATECH ShapeGuard brushes for composite
 (1 x pre-polisher and 1 high lustre)



Restoration Highlight Kit Tips

REF 60020091
 30 x BRILLIANT EverGlow Tips 0.2g, 10 tips of each,
 1 x One Coat 7 Universal Bond 5ml +50 brushes,
 2 x DIATECH ShapeGuard brushes for composite
 (1 x pre-polisher and 1 high lustre)



Buy either kit and get
**A FREE COMPOSITE
 INSTRUMENT SET**

UK promotion only - copy invoice required send to info.uk@coltene.com



Universal submicron hybrid composite

BRILLIANT EverGlow®

Enduring gloss – made brilliant



Thanks to its sophisticated load composition, BRILLIANT EverGlow immediately produces an exceptionally smooth, glossy and satin-finished surface as soon as the filling is put in place. Polishing is simplified and a very aesthetic restoration can be done in no time. In addition, the versatile BRILLIANT EverGlow filling material excels in its enhanced gloss retention.

- Exceptional polishability and long-lasting brilliance
- Aesthetic single-shade restorations
- Versatile shade system with three translucency levels
- Excellent sculptability and smooth consistency
- Good wettability with minimal stickiness to the instrument



BRILLIANT EverGlow® Refills

1x Syringe, 3g OR 20x Tip, 0,2

	BLEACH	Syringe > 60019697 Tips > 60019706
	A1/B1	Syringe > 60019698 Tips > 60019707
	A2/B2	Syringe > 60019699 Tips > 60019708
	A3/D3	Syringe > 60019700 Tips > 60019709
	A3.5/B3	Syringe > 60019701 Tips > 60019710
	C2/C3	Syringe > 60019702 Tips > 60019716
	A4/C4	Syringe > 60019703 Tips > 60019717

	Translucent	Syringe > 60019704 Tips > 60019718
	Bleach Translucent	Syringe > 60019705 Tips > 60019719
	OPAQUE A3	Syringe > 60019745 Tips > 60019748
	OPAQUE A1	Syringe > 60019744 Tips > 60019747
	OPAQUE BLEACH	Syringe > 60019743 Tips > 60019746

WATCH THE VIDEOS



CLICK OR SCAN

BRILLIANT EverGlow®



BRILLIANT EverGlow® System Kit Syringes

REF 60019722

9 x syringes, 3g 1 x BL, A1/B1, A2/B2, A3/D3, A3.5/B3, C2/C3, A4/C4, BL Trans, Translucent & Shade guide

BRILLIANT EverGlow® System Kit Tips

REF 60019723

60 x Tips 0.2 g - 5 x BL, 8 x A1/B1, A2/B2, A3/D3, A3.5/D3, 5 x C2/C3, A4/C4, BL Trans, 8 x Translucent & Shade guide

Buy either kit and get
A FREE COMPOSITE INSTRUMENT SET

UK PROMOTION ONLY

Copy invoice required - send to info.uk@coltene.com



BRILLIANT EverGlow® Flow

BRILLIANT EverGlow Flow provides convincing aesthetic results and is the perfect partner for the universal BRILLIANT EverGlow composite. The aesthetic composite is easy to handle, and suitable for a host of different indications.

BRILLIANT EverGlow® Flow

BRILLIANT EverGlow® Flow

1 x 2g Syringe

5 x Standard application needles

5 x Fine application needles 0.4 mm




	BLEACH BL	Syringe > 60019752
	A1/B1	Syringe > 60019753
	A2/B2	Syringe > 60019754
	A3/D3	Syringe > 60019755
	A3.5/B3	Syringe > 60019756
	A4/C4	Syringe > 60019757
	TRANSLUCENT TRANS	Syringe > 60019758
	OPAQUE A2 OA2	Syringe > 60019751

Reinforced CAD/CAM composite for permanent restoration

BRILLIANT Crios Bloc

BRILLIANT Crios is a high-performance composite bloc for permanent restorations. The material combines all the advantages of an innovative submicron hybrid composite material with those of a CAD/CAM fabrication process for the reliable, aesthetic and fast manufacturing of restorations without a separate firing process.

- High flexural strength for resistant restorations
- Tooth-like modulus of elasticity for a shock-absorbing effect and pleasant bite experience
- Blends in extremely well to give natural aesthetics
- High resistance and low abrasion on opposite tooth
- No glazing by furnace required.
- Can be modified and repaired
- Exceptional grinding precision for more freedom of design
- Easy polishing for fast shine
- Secure bonding due to ONE COAT 7 UNIVERSAL
- The ideal restoration material for every situation



CEREC/inLab – 5 blocs
5 x Blocs (size 12 or 14)

	Low Translucent		High Translucent	
	5 x 12	5 x 14	5 x 12	5 x 14
BLEACH	60019998 BL LT 12	60019985 BL LT 14	-	-
A1	60019999 A1 LT 12	60019986 A1 LT 14	60020007 A1 HT 12	60019994 A1 HT 14
A2	60020000 A2 LT 12	60019987 A2 LT 14	60020008 A2 HT 12	60019995 A2 HT 14
A3	60020001 A3 LT 12	60019988 A3 LT 14	60020009 A3 HT 12	60019996 A3 HT 14
A3.5	60020002 A3.5 LT 12	60019989 A3.5 LT 14	-	-
B1	60020003 B1 LT 12	60019990 B1 LT 14	60020010 B1 HT 12	60019997 B1 HT 14
B2	60020004 B2 LT 12	60019991 B2 LT 14	-	-
B3	60020005 B3 LT 12	60019992 B3 LT 14	-	-
C2	60020006 C2 LT 12	60019993 C2 LT 14	-	-

BRILLIANT Crios Disc

- Made from our award winning Brilliant Everglow composite
- High surface lustre for smoother restorations
- Duo shades for economy in high or low translucencies



The disc is available in the two different heights 14mm/18mm and following shades and translucencies:

	Low Translucent		High Translucent	
	14 mm	18mm	14mm	18mm
A1	60022899 A1 LT H14	60022886 A1 LT H18	6002907 A1 HT H14	60022894 A1 HT H18
A2	60022900 A2 LT H14	60022887 A2 LT H18	60022908 A2 HT H14	60022895 A2 HT H18
A3	60022901 A3 LT H14	60022888 A3 LT H18	60022909 A3 HT H14	60022896 A3 HT H18
A3.5	60022902 A3.5 LT H14	60022889 A3.5 LT H18	-	-

**FREE* One Coat 7
Single Dose Kit**
REF 60019540
with each disc

LEARN MORE



CLICK OR SCAN

Direct composite veneering system

BRILLIANT COMPONEER®



A Smile Made Brilliant –
Reliable. Affordable. Time-saving.

- Aesthetic, high-quality anterior restorations
- Minimally invasive restorations based on composites
- Only one session - no laboratory necessary
- Long-lasting gloss and excellent polishability
- Highly cost-effective for dental practice and patient

LEARN MORE



CLICK OR SCAN

BRILLIANT COMPONEER Intro Kit Syringes


REF 60023214

- 6 × Veneers Upper Jaw, Bleach Translucent, M, 11-13/21-23
- 6 × Veneers Upper Jaw, Bleach Translucent, L, 11-13/21-23
- 6 × Veneers Upper Jaw, Translucent, M, 11-13/21-23
- 6 × Veneers Upper Jaw, Translucent, L, 11-13/21-23
- 2 × Veneers Upper Jaw Premolars, Bleach Translucent, L, 14/24
- 2 × Veneers Upper Jaw Premolars, Translucent, L, 14/24
- 25 × COMPONEER Brushes
- 1 × COMPONEER Contour Guide - Lower Medium
- 1 × Placer
- 1 × COMPONEER Contour Guide - Lower Small
- 1 × Etchant Gel S, 2.5 ml
- 1 × Modelling Instrument MBC
- 2 × Placer Adapter
- 10 × Placer Stamp
- 1 × Silicone polisher 2303RA
- 1 × Silicone polisher 2311RA
- 1 × Silicone polisher 2403RA
- 1 × Silicone polisher 2411RA
- 1 × COMPONEER Contour Guide - Upper Large
- 1 × COMPONEER Contour Guide - Upper Medium
- 1 × COMPONEER Contour Guide - Upper Small
- 1 × BRILLIANT EverGlow A1/B1, 3 g
- 1 × BRILLIANT EverGlow A2/B2, 3 g
- 1 × BRILLIANT EverGlow Bleach Translucent (Trans), 3 g
- 1 × BRILLIANT EverGlow Bleach Translucent (BL Trans), 3 g
- 1 × ONE COAT 7 UNIVERSAL, 5 ml
- 20 × Application needles
- 1 × BRILLIANT COMPONEER Shade Guide

BRILLIANT COMPONEER Intro Kit Tips

REF 60023215

- 6 × Veneers Upper Jaw, Bleach Translucent, M, 11-13/21-23
- 6 × Veneers Upper Jaw, Bleach Translucent, L, 11-13/21-23
- 6 × Veneers Upper Jaw, Translucent, M, 11-13/21-23
- 6 × Veneers Upper Jaw, Translucent, L, 11-13/21-23
- 2 × Veneers Upper Jaw Premolars, Bleach Translucent, L, 14/24
- 2 × Veneers Upper Jaw Premolars, Translucent, L, 14/24
- 25 × COMPONEER Brushes
- 1 × COMPONEER Contour Guide - Lower Medium
- 1 × Placer
- 1 × COMPONEER Contour Guide - Lower Small
- 1 × Etchant Gel S, 2.5 ml
- 1 × Modelling Instrument MBC
- 2 × Placer Adapter
- 10 × Placer Stamp
- 1 × Silicone polisher 2303RA
- 1 × Silicone polisher 2311RA
- 1 × Silicone polisher 2403RA
- 1 × Silicone polisher 2411RA
- 1 × COMPONEER Contour Guide - Upper Large
- 1 × COMPONEER Contour Guide - Upper Medium
- 1 × COMPONEER Contour Guide - Upper Small
- 8 × BRILLIANT EverGlow A1/B1, 0.2 g
- 8 × BRILLIANT EverGlow A2/B2, 0.2 g
- 8 × BRILLIANT EverGlow Translucent (Trans), 0.2 g
- 8 × BRILLIANT EverGlow Bleach Translucent (BL Trans), 0.2 g
- 1 × ONE COAT 7 UNIVERSAL, 5 ml
- 20 × Application needles
- 1 × BRILLIANT COMPONEER Shade Guide



SPECIAL PRICE
INTRO KIT SYRINGES - ONLY £985
INTRO KIT TIPS - ONLY £895

**CONTACT YOUR TERRITORY
MANAGER FOR A DEMONSTRATION**



ONE COAT 7 UNIVERSAL

- Self Etch, Selective Etch and Total Etch, a bond for all adhesive techniques
- With activator for chemically curing products
- Exceptional bonding values on enamel and dentine
- Excellent marginal seal and outstanding marginal leakage tightness
- Ergonomically shaped, triangular bottle with dropper for fine dosing



ONE COAT 7 UNIVERSAL Refill
REF 60019539

Contains MDP - One bond for all indications

LEARN MORE



CLICK OR SCAN

ONE COAT 7 UNIVERSAL Intro Kit

- REF 60019538
- 1 x 5 ml ONE COAT 7 UNIVERSAL
 - 1 x 2.5 ml Etchant Gel S
 - 1 x 20 pcs. Application needles
 - 1 x 50 pcs. Brushes black
 - 1 x Mixing well with cover



UNIVERSAL ACTIVATOR 3ml
for self and dual cure materials
REF 7054
1 x 3ml Activator,
1 x 50pcs. Microblush, black

REQUEST A FREE SAMPLE*

(samples in single dose format)

Self Adhesive Dual Curing Resin Cement



RESTORATION



USER REPORT



CLICK OR SCAN

SoloCem[®]

Secure cementing does not depend on the number of working steps. Self-adhesive SoloCem achieves reliable values - with no additional bonding. This is the result of the sophisticated monomers contained in the composite-like formulation. The convenient automix syringe and small tips make SoloCem a simple and time saving universal cement.

- A true self-adhesive cement: Etching, bonding and cementation in one step
- Long lasting cementation
- Fast and easy handling without requiring a separate adhesive
- Provides high bond strength on a variety of materials
- Unique combination of MDP & 4 - META
- 3 shades: Translucent, Dentine and White Opaque



Enhance your bonding when used together

3+1

BUY 3 - GET 1 FREE*
ON REFILLS

SoloCem Refill Dentin REF 60014062

1 × 5ml / 10.5g SoloCem Dentin
5 × Mixing Tip Brown Short Super Fine
5 × Mixing Tip Brown Short Fine

SoloCem Refill White Opaque REF 60014061

1 × 5ml / 10.5g SoloCem White Opaque / as above

SoloCem Refill Translucent REF 60019840

1 × 5ml / 10.5g SoloCem Translucent / as above

A-silicone impression material

AFFINIS®

The crowning in level of detail



Still one of the fastest A-silicones

AFFINIS fast, with a maximum total setting time of 2 min 10 sec guarantees optimum efficiency and saves time.

Precise moulding of every detail even under the most difficult conditions


The highly active surfactant system is characterized by excellent surface affinity in a moist environment.

Confidence during impression taking

Self-contouring consistency prevents air bubbles and distortions.

Easy monitoring for success

AFFINIS PRECIOUS wash materials with impressive gold and silver pigmentation for excellent readability and detail reproduction.



AFFINIS Putties

includes base 300 ml plus catalyst 300 ml

- REF 6530 putty soft
- REF 6531 fast putty soft
- REF 6535 putty super soft



AFFINIS MonoBody/heavy body

includes 2 x 50 ml cartridges plus 6 x mixing tips

- REF 6590 heavy body
- REF 6650 fast heavy body

includes 2 x 75 ml cartridges plus 8 x mixing tips

- REF 6547 MonoBody
- REF 6520 heavy body
- REF 6620 fast heavy body



AFFINIS light body/regular body

includes 2 x 50 ml cartridges plus 12 x mixing tips

- REF 6501 light body
- REF 6601 fast light body
- REF 6511 regular body
- REF 6611 fast regular body
- REF 6775 PRECIOUS light body
- REF 6776 PRECIOUS regular body

3+1

BUY 3 - GET 1 FREE*
ON ALL REFILLS



AFFINIS® System 360

includes 2 x 380 ml cartridges

- REF 6487 System 360 heavy body
- REF 6482 System 360 Monobody
- REF 6472 Putty



CLICK OR SCAN

WATCH THE VIDEO


BLACK

VISIBLE PRECISION WITH LOVE FOR DETAILS

AFFINIS® BLACK EDITION

The new black tray material AFFINIS heavy body BLACK EDITION reaches with gold or silver coloured AFFINIS PRECIOUS for excellent readability of the impression.

The unique colour combination displays details effectively and precisely which supports the qualitative assessment of the impression result.



AFFINIS BLACK EDITION

includes 2 x 75 ml cartridges plus 8 x mixing tips

REF 60019554 heavy body

includes 2 x 380 ml cartridges

REF 60019776 System 360 heavy body



REF 60019554

“The best impressions the lab had seen”

Dr Mahesh Vasireddy is Principal Dentist at York Dental Practice. He is an experienced dental surgeon, with a special interest in endodontics and cosmetic dentistry. Mahesh believes in providing the highest standard of ethical dentistry and putting his patients first. Here, he discusses why he uses COLTENE's AFFINIS™ Black Edition impression material.


“I've sent my impressions taken with AFFINIS™ to several of my labs - one said they were the best they'd ever seen! The main reason is the black and gold colour combination, so any drag or voids can be picked up straight away, with clear margins.

I started using AFFINIS™ Black Edition heavy body, with AFFINIS™ Precious regular body a couple of years ago. It gives extremely accurate impressions, without distortion and voids. These are fantastic impressions of a full arch prep or a single unit, with excellent colour contrast, that are easy to remove, meaning a comfortable experience for my patients.”



Dr Mahesh Vasireddy
BDS, MFDS (Edin),
MFDS (London), IQE,
MDS (Endodontics, India)

A-Silicone impression material



PRESIDENT The Original

**A SMALL STEP FOR DENTISTS.
A GREAT STEP FOR DENTISTRY.**

- High tear strength and elasticity
- Optimised flowing behavior and thixotropy
- Easier-to-mix putties
- New and more contrasting colours



Putty Putty soft Putty soft fast Putty super soft



Regular body Light body Xtra light body

Putty

- Putty 60019942
 - Putty soft 60019943
 - Putty soft fast 60022166
 - Putty super soft 60019944
- 2 × 300 ml
2 × spoons

Heavy Body

- Heavy body 50 ml 60019940
- Heavy body 75 ml 60019941

Wash cartridge 50 ml

- Regular body 60019939
 - Light body 60019938
 - Xtra light body 60019937
- 2 × 50 ml
12 × mixing tips universal



System 360

- MonoBody 60019946
 - Heavy body 60019948
- 2 × 380 ml



3+1

BUY 3 - GET 1 FREE*
ON ALL REFILLS

Speedex®

The C-silicone with A-silicone quality

- One for all: the Speedex Universal Activator for all wash and tray materials
- New Speedex putty soft material impresses with outstanding final hardness
- You can customise the working time by increasing or reducing the amount of activator paste



3 + 1

BUY 3 REFILLS
GET 1 REFILL FREE*

Speedex Refills

REF 4970
Speedex putty 1 × 910 ml

REF 4980
Speedex light body 1 × 140 ml

REF 4986
Speedex medium 1 × 140 ml

REF 4990
Speedex Universal Activator 1 × 60 ml

Laboratory Modelling Silicone

Lab Putty

Lab-Putty is the outstanding polysiloxane laboratory modelling product material with 63 - 70 Shore A for model repairs, isolation and denture heat protection, keys and block outs and matrixes for formation of general occlusion and functional patterns.

A smooth acrylic surface is guaranteed with the use of Lab Putty hard for relining, keys and repair models. Separate isolation is eliminated. Its high heat-resistance (300° C) and perfect dimensional stability makes it suitable for soldering models (90 Shore A).



3 + 1

BUY 3 TUBS
GET 1 TUB FREE*

Lab Putty
1.3 kg
REF 8805
Single Pack

1 x 900 ml Base
1 x 40 ml Activator
1 x Dosage spoon

Multilayer diamond instruments

DIATECH® Kits



COMPOSITE PREPARATION AND FINISHING KIT

REF 250059AA

This kit comprises 12 Multilayer FG Diamonds in 2 different grit sizes, 2 Silicone polishers and 1 polishing brush for contra-angles. The diamond coated cutting surface of these instruments is pointed at an angle of approximately 95°. The range of applications of these diamonds covers trimming, contouring and finishing the occlusal surfaces of composite restorations and finishing porcelain inlays and onlays chairside. The polisher as well as the polishing brushes provide the surface with a tooth-like enamel gloss.



TOTALLY AESTHETICS KIT By Dr. Monik Vasant

REF 60023603

Now updated to include new products from the Diatech family including DIATECH® ShapeGuard



COMPONEER PREPARATION AND FINISHING KIT by Dr. Mario Besek

REF 60013844

Prep, Finish, Polish! The selected, coarse, blue-coded diamonds are ideally suited for prep-work, the finer red and white coded diamonds for finishing and contouring. The specially designed finest diamond (orange) and the polishers improve the quality of margins and surface gloss.



TOP 13 UK KIT

REF 250052AA
















13 of the most popular bur patterns purchased by UK dentists for Crown and Bridge preparation, Composite contouring and Finishing.

Kit contains 26 burs (2 of each)

DIATECH® Diamond Instruments

TOP SELLERS

BUY ANY 3 PACKS GET 1 PACK FREE!*

		GRIT	PACK SIZE	REF CODE
	G368-314-020-5	F	5 pcs	200010AA
	G368-314-023-5	F	5 pcs	200014AA
	G379-314-023-5	F	5 pcs	200037AA
	G379-314-023-5	XF	5 pcs	200041AA
	G801-314-010	ML	5 pcs	200079AA
	G801-314-012	ML	5 pcs	200081AA
	G801-314-014	ML	5 pcs	200083AA
	G801-314-016	ML	5 pcs	200086AA
	G801-314-018	ML	5 pcs	200089AA
	G835-314-010-4	ML	5 pcs	200203AA
	G835-314-012-4	ML	5 pcs	200206AA
	G859-314-010-10	F	5 pcs	200502AA
	G863-314-012-10	XF	5 pcs	200562AA
	G889-314-009-3.5	F	5 pcs	200732AA
	G889-314-009-3.5	ML	5 pcs	200733AA

Infection Control for every practice

SciCan
OPTIM[®] 1

1 MINUTE
CONTACT TIME**

ONE-STEP CLEANER
& DISINFECTANT

Active Ingredient: 0.5% Hydrogen Peroxide



With OPTIM 1 there is no need to compromise between efficacy and being protected. With a fast broad-spectrum contact time it is one of the fastest and most effective cleaner & disinfectants. It is also safe and gentle for you and surfaces. OPTIM 1 gives you confidence that you are compliant and protecting your patients, staff and the environment; the perfect balance.

Demonstrated virucidal efficacy also extends to new and emerging pathogens, with a **30 second contact time applying to Novel Coronavirus (COVID-19)**, Severe Acute Respiratory Syndrome (SARS CoV) and Middle East Respiratory Syndrome (MERS CoV).



FAST AND GENTLE

1 × FREE DISPENSER PUMP*

One per practice while stocks last. Case sale only. Send copy invoice to info.uk@coltene.com



EFFECTIVE
AGAINST COVID-19

** Effective within 1 minute according to European Standards; EN 14476, EN 14348, EN 13727, EN 13624, EN 16615. According to the required RKI/DVV standard in Germany SciCan OPTIM 1 is FULL VIRUCIDAL within 2 minutes

ACTIVE AGAINST COVID-19



Fast Broad-Spectrum Disinfectant

Surfaces must remain wet for the entire contact time to achieve disinfection.

- ▶ Adeno, Polio, Murino Noro in only **30 seconds** according to EN 14476
- ▶ Tuberculocidal & Fungicidal in **1 minute**
- ▶ Bactericidal & Virucidal in **30 seconds** according to EN 14476
- ▶ **Full virucide in 2 minutes** according to RKI/DVV standard (Polio, Adeno, Murino Noro, SV40)
- ▶ Effective in **1 minute** according to EN 16615
- ▶ Excellent wettability – stays wet long enough for contact times to elapse



Effective One-Step Cleaner

An effective 1 step cleaner disinfectant removes the organic debris while it disinfects.



Compatible

Compatible with a wide range of hard, non-porous materials and surfaces.



Safe and Gentle

- ▶ Free of artificial scents
- ▶ Non-Toxic
- ▶ Non-irritating to skin, eyes and respiratory system
- ▶ No GHS hazard pictograms, signal words or statements
- ▶ Leaves no active residue; doesn't contribute to antimicrobial resistance



Environmentally Friendly

The active ingredient Hydrogen Peroxide breaks down into water and oxygen leaving no active residues which reduces the environmental impact.

SciCan OPT/M 1 Ready-to-use Pre-saturated Wipes

RTU Wipes 15 x 18 cm (12 tubs x 160 sheets per tub) REF OPTEU1-12W
 RTU Wipes 25 x 25 cm (12 tubs x 60 sheets per tub) REF OPTEU1-W10x10

SciCan OPT/M 1 Ready-to-use Liquid

1 case of 4 x 4 Litres REF OPTEU1-4x4
 Pump for 4 Litre bottle REF FND30

**Ask your territory manager for your
 FREE SAMPLE tub of wipes - to see how
 effective the switch to OPT/M 1 can be.**

FOR MORE
 INFORMATION



CLICK OR SCAN

MEET THE TEAM



Debi Boynton
07764 327 425
debi.boynton@coltene.com
North West
Scotland - Postcodes:
DG, FK, G1-G8, K8, ML, PA, PH



Sue Robins
07786 965 019
sue.robins@coltene.com
East Midlands



Chris Nolan
07889 227 424
chris.nolan@coltene.com
North East
Scotland - Postcodes:
AB, DD, EH, HS, IV, KW, KY TD, ZE



Nick O'Keeffe
00353 8641 35766
nicholas.okeeffe@coltene.com
Ireland



Emma Gilks
07906 749 682
emma.gilks@coltene.com
West Midlands & North Wales



Shah Monnan
07712 407 071
shah.monnan@coltene.com
Essex & East London



Andrew Miller
07788 583 057
andrew.miller@coltene.com
South Wales & West Country



Kelly McCormack
07818 424 210
kelly.mccormack@coltene.com
North, West & Central London



Louise Fish
07788 146 130
louise.fish@coltene.com
South London & South Coast

UK OFFICE - SUSSEX



Mark Allen
07787 570 517
Managing Director



Nic Coomber
01444 235 486 - Ext 226
National Account Manager
& Marketing Strategist



Helen Chart
01444 235 486 - Ext 227
Finance & Office Manager



Alicia Franks
01444 235 486 - Ext 223
Office Coordinator

All promotional offers valid to 30.06.2021 except where stated. *Free stock supplied by dealer unless stated otherwise.

© 2021 COLTENE Printing errors mistakes and changes reserved.